



CHROMATOGRAPHIC PROCESS MEDIA CATALOG



TOSOH BIOSCIENCE

ABOUT US

WITH A GLOBAL PERSPECTIVE.

TOSOH BIOSCIENCE GmbH, Separations Business Unit, Stuttgart, is an acknowledged global leader in the field of bioseparations. Established as TosoHaas in 1987, the original joint venture between Tosoh Corporation of Japan and the Rohm and Haas Company, USA, has become synonymous with advanced products and quality support. In the year 2000, Tosoh Corporation acquired a 100% controlling interest changing the name to TOSOH BIOSEP. In the course of unifying all Tosoh affiliates, the new Brand Name Tosoh Bioscience evolved. Today, the two branches, Bioseparations and Diagnostics operate with the same name Tosoh Bioscience - Separations Business Unit and accordingly Diagnostics Business Unit. Tosoh manufacturing sites in Japan provide products to the sales and support subsidiaries in the U.S. and Europe, ensuring full global coverage. Over the last 30 years, TSKgel SW columns have become the worldwide industry standard for size exclusion chromatography of biomolecules.



1
TOSOH BIOSCIENCE GMBH

ZETTACHRING 6
70567 STUTTGART
GERMANY

T · + 49 (0) 711 · 13257 - 0
F · + 49 (0) 711 · 13257 - 89
INFO.TBG@TOSOH.COM
WWW.TOSOHBIOSCIENCE.COM

2
TOSOH BIOSCIENCE LLC

3604 HORIZON DRIVE,
SUITE 100
KING OF PRUSSIA, PA 19406, USA

T · +1 484-805-1219
F · +1 610-272-3028
INFO.TBL@TOSOH.COM
WWW.TOSOHBIOSCIENCE.COM

3
TOSOH CORPORATION

3-8-2 SHIBA, MINATO-KU
TOKYO 105-8623
JAPAN

T · +81 3-5427-5118
F · +81 3-5427-5198
INFO@TOSOH.CO.JP
WWW.TOSOHBIOSCIENCE.COM



4
TOSOH BIOSCIENCE SHANGHAI CO. LTD

ROOM 301, PLAZA B,
NO. 1289 YI SHAN ROAD
XU HUI DISTRICT
SHANGHAI, 200233, CHINA

T · +86 21-3461-0856
F · +86 21-3461-0858
INFO@TOSOH.COM.CN
WWW.SEPARATIONS.ASIA.TOSOHBIOSCIENCE.COM

5
TOSOH ASIA PTE. LTD

63 MARKET STREET #10-03
SINGAPORE 048942,
SINGAPORE

T · +65 6226-5106
F · +65 6226-5215
INFO@TOSOH.COM.SG
WWW.SEPARATIONS.ASIA.TOSOHBIOSCIENCE.COM

TOSOH HISTORY

- 1935 FOUNDING OF TOYO SODA MANUFACTURING CO., LTD.
- 1936 OPERATION OF NANYO MANUFACTURING COMPLEX BEGINS
- 1971 SCIENTIFIC INSTRUMENTS DIVISION FORMED, FIRST GPC COLUMN USING TSK-GEL DEVELOPED BY TOSOH
- 1974 HIGH PERFORMANCE LIQUID CHROMATOGRAPHY COLUMN PLANT IS COMPLETED
- 1979 TOSOH DEVELOPS TOYOPEARL MEDIA
- 1983 TOSOH DEVELOPS HYDROPHOBIC INTERACTION MEDIA
- 1987 TOSOH US OPERATIONS FORMED IN MONTGOMERYVILLE
- 1989 TOSOH GMBH OPERATIONS FORMED IN STUTTGART
- 1995 TOSOH NANYO GEL FACILITY RECEIVES ISO 9001
- 2000 IN NOVEMBER FORMER TOSOH US OPERATIONS BECOMES TOSOH BIOSEP LLC, A 100% SUBSIDIARY OF TOSOH CORPORATION
- 2001 IN JANUARY FORMER TOSOH GMBH EUROPEAN OPERATIONS BECOMES TOSOH BIOSEP GMBH, A 100% SUBSIDIARY OF TOSOH CORPORATION
- 2002/ TOSOH CORPORATION ANNOUNCES THAT ALL TOSOH AFFILIATED SCIENTIFIC AND DIAGNOSTIC SYSTEM
- 2003 RELATED COMPANIES IN EUROPE, WILL BE UNIFIED UNDER THE NEW NAME TOSOH BIOSCIENCE.
- 2008 ECOSEC , THE 7TH GENERATION GPC SYSTEM IS INTRODUCED GLOBALLY
- 2009 TOSOH BIOSCIENCE GMBH CELEBRATES ITS 20TH ANNIVERSARY IN STUTTGART
- 2010 TOSOH CELEBRATES ITS 75TH YEAR IN BUSINESS WITH THE OPENING OF FIVE NEW PLANTS, AND CONTINUED RAPID EXPANSION IN CHINA
- 2011 AFTER DEVELOPING THE FIRST TSK_{gel} GPC COLUMN IN 1971, TOSOH NOW LOOKS BACK ON 40 YEARS EXPERIENCE IN SUCCESSFUL TECHNOLOGY IN SIZE EXCLUSION AND GEL PERMEATION CHROMATOGRAPHY

1
TOSOH BIOSCIENCE GMBH

ZETTACHRING 6
70567 STUTTGART
GERMANY

T · + 49 (0) 711 · 13257 - 0
F · + 49 (0) 711 · 13257 - 89
INFO.TBG@TOSOH.COM
WWW.TOSOHBIOSCIENCE.COM

2
TOSOH BIOSCIENCE LLC

3604 HORIZON DRIVE,
SUITE 100
KING OF PRUSSIA, PA 19406, USA

T · +1 484-805-1219
F · +1 610-272-3028
INFO.TBL@TOSOH.COM
WWW.TOSOHBIOSCIENCE.COM

3
TOSOH CORPORATION

3-8-2 SHIBA, MINATO-KU
TOKYO 105-8623
JAPAN

T · +81 3-5427-5118
F · +81 3-5427-5198
INFO@TOSOH.CO.JP
WWW.TOSOHBIOSCIENCE.COM



4
TOSOH BIOSCIENCE SHANGHAI CO. LTD

ROOM 301, PLAZA B,
NO. 1289 YI SHAN ROAD
XU HUI DISTRICT
SHANGHAI, 200233, CHINA

T · +86 21-3461-0856
F · +86 21-3461-0858
INFO@TOSOH.COM.CN
WWW.SEPARATIONS.ASIA.TOSOHBIOSCIENCE.COM

5
TOSOH ASIA PTE. LTD

63 MARKET STREET #10-03
SINGAPORE 048942,
SINGAPORE

T · +65 6226-5106
F · +65 6226-5215
INFO@TOSOH.COM.SG
WWW.SEPARATIONS.ASIA.TOSOHBIOSCIENCE.COM

TOSOH HISTORY

- **1935** FOUNDING OF TOYO SODA MANUFACTURING CO., LTD.
- 1936** OPERATION OF NANYO MANUFACTURING COMPLEX BEGINS
- 1971** SCIENTIFIC INSTRUMENTS DIVISION FORMED, FIRST GPC COLUMN USING TSK-GEL DEVELOPED BY TOSOH
- 1974** HIGH PERFORMANCE LIQUID CHROMATOGRAPHY COLUMN PLANT IS COMPLETED
- 1979** TOSOH DEVELOPS TOYOPEARL MEDIA
- 1983** TOSOH DEVELOPS HYDROPHOBIC INTERACTION MEDIA
- 1987** TOSOH US OPERATIONS FORMED IN MONTGOMERYVILLE
- 1989** TOSOH GMBH OPERATIONS FORMED IN STUTTGART
- 1995** TOSOH NANYO GEL FACILITY RECEIVES ISO 9001
- 2000** IN NOVEMBER FORMER TOSOH US OPERATIONS BECOMES TOSOH BIOSEP LLC, A 100% SUBSIDIARY OF TOSOH CORPORATION
- 2001** IN JANUARY FORMER TOSOH GMBH EUROPEAN OPERATIONS BECOMES TOSOH BIOSEP GMBH, A 100% SUBSIDIARY OF TOSOH CORPORATION
- 2002/** TOSOH CORPORATION ANNOUNCES THAT ALL TOSOH AFFILIATED SCIENTIFIC AND DIAGNOSTIC SYSTEM
- 2003** RELATED COMPANIES IN EUROPE, WILL BE UNIFIED UNDER THE NEW NAME TOSOH BIOSCIENCE.
- 2008** ECOSEC , THE 7TH GENERATION GPC SYSTEM IS INTRODUCED GLOBALLY
- 2009** TOSOH BIOSCIENCE GMBH CELEBRATES ITS 20TH ANNIVERSARY IN STUTTGART
- 2010** TOSOH CELEBRATES ITS 75TH YEAR IN BUSINESS WITH THE OPENING OF FIVE NEW PLANTS, AND CONTINUED RAPID EXPANSION IN CHINA
- 2011** AFTER DEVELOPING THE FIRST TSKgel GPC COLUMN IN 1971, TOSOH NOW LOOKS BACK ON 40 YEARS EXPERIENCE IN SUCCESSFUL TECHNOLOGY IN SIZE EXCLUSION AND GEL PERMEATION CHROMATOGRAPHY



CONTENTS



➤		ABOUT US
	page [00 - 01]	CONTENTS
	page [02 - 06]	GENERAL INTRODUCTION
	page [07 - 07]	WHAT'S NEW

➤	SEC	SIZE EXCLUSION CHROMATOGRAPHY
	page [08 - 14]	ORDERING INFORMATION
	page [15 - 15]	

➤	IEC	ION EXCHANGE CHROMATOGRAPHY
	page [16 - 25]	ORDERING INFORMATION
	page [26 - 29]	

➤	HIC	HYDROPHOBIC INTERACTION CHROMATOGRAPHY
	page [30 - 39]	ORDERING INFORMATION
	page [40 - 41]	

➤	AFC	AFFINITY CHROMATOGRAPHY
	page [42 - 49]	ORDERING INFORMATION
	page [50 - 51]	

➤	CUSTOM DESIGNED RESINS	
	page [52 - 56]	CUSTOM DESIGNED RESINS



INTRODUCTION

OUR HISTORY

Tosoh Bioscience is an acknowledged leader in the field of bioseparations. Established in August 1987 as TosoHaas, we were the joint venture of two multi-billion dollar companies: Tosoh Corporation of Japan and Rohm and Haas Company, USA. In the year 2000, Tosoh Corporation purchased Rohm and Haas Company's business interest in TosoHaas and formed Tosoh Biosep subsidiaries with sales, marketing, and technical service offices in the U.S. and Germany. In 2002, Tosoh Biosep was re-organized with Tosoh's clinical diagnostic business units under the name of Tosoh Bioscience. In 2007, Tosoh Scientific Instrument Division changed its name to Tosoh Bioscience Division and Tosoh (Shanghai) Co., Ltd. was also created.

A GLOBAL PERSPECTIVE

Since our beginnings, we have become an effective manufacturing partner with many of the world's most innovative pharmaceutical and biotechnology companies. Our TSKgel® and Toyopearl® chromatographic resins are important components in the analysis, isolation, and purification of biomolecules. Tosoh's portfolio of over 500 specialty items encompassing all common modes of liquid chromatography can help purify any protein, peptide, enzyme, nucleic acid, antibiotic, or small molecule that may be encountered.

SEAMLESS SCALE-UP WITH POLYMER RESINS

The polymeric TSKgel PW-type analytical HPLC columns and Toyopearl production scale resins are made using the same chemistry. This gives them the same chromatographic selectivity so that an analytical method developed on a TSKgel column can be easily scaled to the corresponding bulk Toyopearl resin for manufacturing use. This product synergy saves valuable development time.

TOSOH'S TECHNOLOGY

For over thirty-five years our parent, Tosoh Corporation, has been a world leader in the analysis and purification of proteins. Tosoh's TSKgel SW columns are the industry standard for size exclusion chromatography of biomolecules. This has led Tosoh to a fundamental understanding of the role played by pore diameter and molecular size in chromatographic separations. This knowledge allows Tosoh to design higher performance polymeric resins for ion exchange, hydrophobic interaction, affinity and reversed phase applications.

NEW PRODUCT CREATION

Tosoh Bioscience and Tosoh provide solutions for today's biological purification needs. The Custom Resin section of this catalog describes how we are ready to solve future separation challenges. Whether it is a surface modification of an existing product or the creation of a new resin, we encourage you to contact us to provide you with a customized solution.



TOSOH BIOSCIENCE GMBH ZETTACHRING 6 70567 STUTTGART
GERMANY T · + 49 (0) 711 · 13257 - 0 F · + 49 (0) 711 · 13257 - 89
INFO.TBG@TOSOH.COM WWW.TOSOHBIOSCIENCE.COM



TOSOH BIOSCIENCE LLC 3604 HORIZON DRIVE, SUITE 100, KING OF
PRUSSIA, PA 19406, USA T · +1 484-805-1219, F · +1 610-272-3028
INFO.TBL@TOSOH.COM WWW.TOSOHBIOSCIENCE.COM



TOSOH CORPORATION - NANYO COMPLEX
TOYOPEARL RESINS ARE PRODUCED AT THE TOSOH NANYO
COMPLEX – A SPRAWLING 3 MILLION SQ. METER FACILITY THAT IS
JAPAN'S LARGEST CHEMICAL MANUFACTURING COMPLEX.

TOSOH BIOSCIENCE SHANGHAI CO., LTD.
ROOM 301, PLAZA B, NO. 1289 YI SHAN ROAD, XU HUI DISTRICT
SHANGHAI, 200233, CHINA, T · +86 21-3461-0856 F · +86 21-3461-0858
INFO@TOSOH.COM.CN WWW.TOSOHBIOSCIENCE.COM

TOSOH ASIA PTE. LTD
63 MARKET STREET #10-03, SINGAPORE 048942, SINGAPORE,
T · +65 6226-5106 F · +65 6226-5215 INFO@TOSOH.COM.SG
WWW.TOSOHBIOSCIENCE.COM



INTRODUCTION

THE SKILLS TOOLBOX

The science of biopurification is one of the fastest growing technologies in the world. Companies that accept the pace and are flexible will win. At Tosoh, we engineer: particle size, pore size, surface area, ligand chemistry, and innovation. The “Tosoh Resin Innovation Program” (TRIP) discussed in the “Custom Resin” section of this catalog describes our core technologies for solving customer problems.

CUSTOMER TRUST IS OUR NO.1 PRIORITY

Customer trust is earned, not given. Since 1987, Tosoh Bioscience has helped develop processes from benchtop discovery to multi-thousand liter manufacturing scale. A number of our process customers have been with us continuously for over 20 years. We have earned their trust. From development to new product launch we share in the excitement of our customers’ successes, as well as their clinical disappointments. In either case, Tosoh Bioscience uses its technology, innovation, and customer commitment to expedite the development process to build better customer relationships and trust for the future.

Let’s start this process today!

WE LISTEN TO EVERY CUSTOMER REQUEST

Developing trust with our customers allows us to have very detailed discussions about individual customer needs. These discussions guide us in our understanding of emerging new product requirements and also the level of satisfaction the customer has with us as a company. We specialize in designing new products based on existing chemistries to solve specific customer applications.

Do you need a better HETP?
We can do it.

More dynamic binding capacity?
We can do it.

A new ligand immobilization?
We can do it.

Just ask us.
We will listen.

THE BEST EXPERIENCE FOR OUR CUSTOMERS

As your manufacturing partner we want our logistical support to be seamless to you. From your initial product evaluation to full scale manufacturing, our technical support, order administration, and final shipment procedures are in place to make your experience with us a thoroughly enjoyable one. Courteous professionals are available to answer all of your questions. Our customer service processes are also in place to deliver our resins to you as agreed. Journal references, product specifications, samples, quotations, applications and scale up support are all just a telephone call, email, or mouse-click away.

Tosoh Bioscience works hard to make it easy for you to do business with us.

COMPLETE REGULATORY SUPPORT

Controlled production processes and auditable manufacturing facilities are the regulatory requirements of our industry. Well characterized, consistent resin products are as important to us as your therapeutic molecules with their known impurity profiles are to you. For this reason Tosoh Corporation’s manufacturing site in Japan and both Tosoh Bioscience subsidiaries in the US and Europe are ISO 9001:2000 certified. Our manufacturing facilities are continuously audited by customers and we pass their tests. We maintain regulatory support information on our chromatographic resins available to our customers for their process registration needs.

Our resin processes are well controlled.

COMPREHENSIVE TECHNICAL SUPPORT

Do you need help developing an HPLC assay for the analysis of a new therapeutic target? How about monitoring the therapeutic’s metabolites in the human body? Do you need help with a “self-packing” process column or a multi-thousand liter one? Our technical support specialists can provide assistance in all of these areas and more. We offer on-site training, application specific seminars, and are the sole sponsor of the HIC/RPC Bioseparation Conference series (www.HIC-RPC.org).

Do you need some help or training? Call us.
We’ll be there for you.



INTRODUCTION

THE PRODUCTS TOOLBOX

From the research laboratory to full scale manufacturing, we offer the same polymer chemistries in our TSKgel and Toyopearl products (Table I). Whether you are scaling up from a TSKgel column HPLC method to Toyopearl for manufacturing, or are scaling down from a Toyopearl purification to the corresponding TSKgel column for the QC of your target, we make it easy to develop methods to do both.

TSKgel COLUMNS

Our TSKgel prepacked columns for high performance liquid chromatography are used for the analysis and purification of proteins, peptides, biopolymers and low molecular weight compounds. Many different chromatographic modes are offered including: size exclusion, ion exchange, hydrophobic interaction, reversed phase, affinity and normal phase chromatography. The column packing materials are either spherical silica or polymer particles, ranging in size from 2 μm to 30 μm . Columns are constructed of stainless steel, glass, or polyetheretherketone (PEEK) and are available from microbore to preparative sizes. Please see our website for more information.

TSKgel RESINS

The polymeric resins with particle sizes of 20 μm and 30 μm used in TSKgel columns are also available in bulk quantities for large scale ion exchange and hydrophobic interaction chromatography. Their mechanical stability and permeability make them excellent for use when increased separation performance and plate count are needed for optimum preparative or process chromatography.

TOYOPEARL RESINS

Toyopearl resins are hydrophilic macroporous methacrylic resins for large-scale chromatographic applications. Their rigid polymeric backbone has better pressure-flow properties than most other commercially made materials. Therefore, higher linear operating velocities can be used for faster process throughput and decreased recycling times. Toyopearl resins are stable over the pH 2-12 range for normal operating conditions and pH 1-13 for cleaning conditions. The resins are available in average particle sizes of 35 μm , 65 μm , 75 μm , and 100 μm for high resolution, intermediate purification, or capture chromatography. Toyopearl resins are also offered in many different pore diameters for size exclusion, ion exchange, hydrophobic interaction, and affinity chromatography.

TABLE I

COMPARISON OF PREPACKED COLUMNS WITH RESPECTIVE BULK POLYMERIC MEDIA

MODE	TSKgel COLUMN FOR ANALYSIS	TOYOSCREEN PROCESS DEVELOPMENT COLUMNS	BULK POLYMERIC MEDIA FOR SCALE-UP AND PRODUCTION	
SEC	TSKgel G1000PW - G6000PW		Toyopearl HW-40 - HW-75	
IEC	TSKgel SP-5PW	ToyoScreen SP-650M	TSKgel SP-5PW	Toyopearl SP-650S, M or C
		ToyoScreen SP-550C		Toyopearl SP-550C
		ToyoScreen MegaCap II SP-550EC		Toyopearl MegaCap II SP-550EC
		ToyoScreen GigaCap S-650M		Toyopearl GigaCap S-650M
		ToyoScreen GigaCap CM-650M		Toyopearl GigaCap CM-650M
	TSKgel CM-5PW	ToyoScreen CM-650M		Toyopearl CM-650S, M or C
	TSKgel DEAE-5PW	ToyoScreen DEAE-650M	TSKgel DEAE-5PW	Toyopearl DEAE-650S, M or C
	TSKgel SuperQ-5PW	ToyoScreen SuperQ-650M	TSKgel SuperQ-5PW	Toyopearl SuperQ-650S, M or C
		ToyoScreen GigaCap Q-650M		Toyopearl GigaCap Q-650M
		ToyoScreen QAE-550C		Toyopearl QAE-550C
ToyoScreen Q-600C AR			Toyopearl Q-600C AR	
HIC	TSKgel Ether-5PW	ToyoScreen Ether-650M	TSKgel Ether-5PW	Toyopearl Ether-650S, M
	TSKgel Phenyl-5PW	ToyoScreen Phenyl-650M	TSKgel Phenyl-5PW	Toyopearl Phenyl-650S, M or C
		ToyoScreen Phenyl-600M		Toyopearl Phenyl-600M
		ToyoScreen Butyl-650M		Toyopearl Butyl-650S, M or C
		ToyoScreen Butyl-600M		Toyopearl Butyl-600M
		ToyoScreen Hexyl-650C		Toyopearl Hexyl-650C
		ToyoScreen PPG-600M		Toyopearl PPG-600M
	ToyoScreen SuperButyl-550C		Toyopearl SuperButyl-550C	
AFC	TSKgel Chelate-5PW	ToyoScreen AF-Chelate-650M		Toyopearl AF-Chelate-650M
			TSKgel Tresyl-5PW	Toyopearl AF-Tresyl-650M
		ToyoScreen AF-Blue HC-650M		Toyopearl AF-Blue HC-650M
		ToyoScreen AF-Red-650M		Toyopearl AF-Red-650M
		ToyoScreen AF-HeparinHC-650M		Toyopearl AF-HeparinHC-650M
	ToyoScreen AF-rProtein-650F		Toyopearl AF-rProtein-650F	

INTRODUCTION



Pore diameter and surface area can be optimized to ensure excellent kinetic access and binding capacity of your target regardless of molecular size. For predictable results in scale-up, Toyopearl resins are based on the same chemistries as the prepacked TSKgel columns. This allows the seamless direct scale-up of methods developed on TSKgel products to Toyopearl.

TOYOSCREEN PROCESS DEVELOPMENT COLUMNS

ToyoScreen process development products are easy to use. They are available as 1 mL and 5 mL prepacked columns of your favorite Toyopearl hydrophobic interaction, ion exchange or affinity resins. They can be connected to most laboratory chromatographic systems for early development resin screening.





INTRODUCTION

SAFETY DATA AND WARRANTY

Tosoh Bioscience maintains Material Safety Data Sheets (MSDS for the U.S. and ERSDS for Europe) on all of its bulk resins. These data sheets contain pertinent information that you may need to protect your employees and customers against any known health or safety hazards associated with our products. Tosoh Bioscience supplies copies of our Material Safety Data Sheets with all bulk resins. You and your customers are responsible for knowing all information and precautions disclosed in the MSDS and any other available materials provided by Tosoh Bioscience which set forth information concerning our products. These describe precautions to be taken in the storage and handling of our products and in the maintenance of the health and safety of persons exposed to our products, the public and the environment with respect to our products. You and your customers should convey such information and precautions to the persons who may be exposed to our products.

We also suggest you contact your supplier of other materials recommended for use with our products for appropriate health and safety precautions prior to their use.

Many of the Toyopearl and TSKgel products have been assigned drug master file numbers by the FDA. Please direct inquiries to the Technical Service Department, Tosoh Bioscience GmbH, Zettachring 6, 70567 Stuttgart, Germany.

Tosoh Bioscience warrants that at the time of delivery each Product conforms to the specifications thereof contained in the Certificate of Analysis ("COA") or the Operating Conditions and Specifications ("OCS") sheet, as relevant, as is provided together therewith; provided, however, that the foregoing warranty applies only if the Products have been properly handled, stored and used by Buyer. Buyer acknowledges that the conditions of use and application by Buyer or its customers of any Product and technical information provided in connection therewith (whether verbal, written or by product evaluation) are beyond Tosoh Bioscience's control. Buyer agrees that it shall test such Product and technical information to determine to its satisfaction whether they are suitable for the intended use and applications of Buyer or its customers. All technical information shall be given without warranty or representation. Buyer assumes, and releases Tosoh Bioscience from, all liabilities, in contract, tort or otherwise, incurred in connection with the use or application of any Product or technical information provided in connection therewith (whether verbal, written or by product evaluation). THIS WARRANTY AND THE PATENT WARRANTY CONTAINED IN THE TERMS AND CONDITIONS ARE GIVEN AND ACCEPTED IN LIEU OF ANY OTHER WARRANTIES AND REPRESENTATIONS, EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION THOSE OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, OR NONINFRINGEMENT OF INTELLECTUAL PROPERTY RIGHTS. Please also see our terms and conditions for other conditions applicable to our warranty.

TECHNICAL DATA AND TRADEMARKS

The technical information and data herein contained (the "Technical Data") are based on information Tosoh Bioscience believes to be reliable and are offered in good faith, but are given without warranty or representation, as the conditions of use and application by you or your customers of our products and the Technical Data are beyond the control of Tosoh Bioscience. You should test our products and Technical Data to determine to your satisfaction whether they will be suitable for the intended use and applications of you or your customers. Suggestions for the uses of our products should not be understood as recommending the use of our products in violation of any patent or other intellectual property right or as permission or license to use any patent or other intellectual property right.

Toyopearl, ToyoScreen, TSK-GEL, TSKgel, Toyopearl MegaCap, and Toyopearl GigaCap are registered trademarks of Tosoh Corporation, the parent company of Tosoh Bioscience.

Triton is a registered trademark of Union Carbide Chemicals and Plastics Co., Inc.

Capto and Source are registered trademarks of GE Healthcare Bio-Sciences.

RETURN POLICY

Tosoh Bioscience is fully committed to delivering quality products and services. TSKgel column products are accompanied by a chromatogram demonstrating the performance of a test mixture on that column and by an OCS sheet that contains information about the Operating Conditions and Specifications of the column. TSKgel and Toyopearl media products are accompanied by Certificates of Analysis.

Despite our commitment to product quality, columns and resins occasionally perform differently than expected in a customer's application. Therefore, we ask you to inspect your TSKgel columns and Toyopearl media within 30 days of receipt by using the conditions described in the documentation supplied with the product. Let us know within this 30 day period if the product does not meet the specifications on the OCS (Operating Conditions and Specifications) sheet and QC document, or as listed in the Certificate of Analysis. Subject to prior authorization, we will accept for return all products that do not perform according to their specifications.

If laboratory TSKgel columns or media containers up to 1 L in size are authorized for return for reasons other than Tosoh Bioscience's error or because of a product defect, there will be a restocking charge of \$ 50 or 10% of the list price, whichever is greater. Most 5 L and 50 L containers of bulk resin are packaged and made to order in Japan and are not returnable except for reasons of Tosoh Bioscience error or of product defect.



WHAT'S NEW

WHAT'S NEW

TOYOPEARL AF-rPROTEIN A-650F

The new Toyopearl AF-rProtein A-650F affinity resin is well-suited for high capacity capturing of immunoglobulins out of high titre feedstocks. It binds human and mouse immunoglobulin G (IgG), as well as Fab fragments. The resin combines a high h-IgG dynamic binding capacity at short residence time and fast mass transfer kinetics with an improved alkaline stability. It provides a typical dynamic binding capacity for IgG of more than 30 mg/ml resin at two minutes residence time. Thus, Toyopearl AF-rProtein A-650F offers the opportunity to improve process economics by increasing both, throughput and robustness.

A newly developed recombinant ligand has been engineered to maintain capacity even after repeated exposure to alkaline solution. Its multipoint attachment to the Toyopearl matrix further enhances chemical and thermal stability. In practice this pays off for a low level of protein A leaching and also for a high resistance to alkaline solutions employed in cleaning-in-place (CIP) procedures.

TOYOPEARL Q-600C AR

Toyopearl Q-600C AR is a new ion exchange resin with an average particle size of 100 µm. It is well-suited for high capacity capture of biotherapeutics or plasma proteins. 'AR' stands for alkaline resistance as the resin combines high dynamic binding capacity and fast mass transfer kinetics with an improved alkaline stability. Thus, Toyopearl Q-600C AR is enabling both, high process throughput and improved process robustness.


Toyopearl Q-600C AR has a slightly larger pore diameter than the well-established Toyopearl QAE-550C resin, but its ligand attachment chemistry creates significantly more dynamic binding capacities of > 100 mg BSA and 90 mg IgG per mL resin. It tolerates caustic cleaning solutions up to 1.0 mol/L NaOH well with essentially no change in dynamic binding capacity (DBC).

TSKgel SP-3PW



The new TSKgel SP-3PW (30) cation exchange resin was designed as a polishing resin with improved binding capacities for peptides and small proteins. It is a strong cation exchange resin having a smaller pore size than the corresponding TSKgel SP-5PW material and also a different selectivity.

Due to an average diameter of 30 µm and a very narrow particle size distribution it is ideally suited when resolution is an issue. Figure 1 shows a microscopic image of TSKgel SP-3PW (30) particles. The resin was developed to provide high DBCs for peptides and small proteins. A typical DBC of about 50 g/L for insulin makes this resin attractive for all peptide purification tasks that involve a cation exchange step.

FIGURE 1 MICROSCOPIC IMAGE OF TSKgel SP-3PW



SIZE EXCLUSION
CHROMATOGRAPHY



SEC SIZE EXCLUSION CHROMATOGRAPHY

SEC PRODUCTS

➤ Toyopearl Resins for SEC

Toyopearl HW-40

Toyopearl HW-50

Toyopearl HW-55


Toyopearl HW-65

Toyopearl HW-75

≡ TOSOH FACT

Tosoh has a long history in size exclusion chromatography (SEC). In 1978 Tosoh first introduced porous silica-based SW columns for the isolation of proteins using LC. These first gels were quickly adopted and referred to as the standard for analytical SEC on FPLC and HPLC systems.

The first Toyopearl process resins for size exclusion chromatography were developed in 1979. Toyopearl HW resins are based on polymethacrylate and are available in various particle and pore sizes. They are used as starting material for the production of all other functionalized Toyopearl resins.





SIZE EXCLUSION CHROMATOGRAPHY

THE ROLE OF SIZE EXCLUSION CHROMATOGRAPHY IN PROCESS PURIFICATION

Size exclusion chromatography (SEC), also known as gel filtration, separates molecules in an aqueous mobile phase according to their physical size in solution as they pass through a porous structure. Molecules with a diameter greater than the largest pores within the resin material are unable to enter the particle. Because they are excluded from the pores they travel quickly through the column and elute first. Smaller molecules, which are able to access pores within the resin particles, permeate a larger accessible volume within the column and are eluted later, in order of decreasing molecular weight.

Because SEC has no adsorption capacity and its separation mechanism dilutes the sample upon elution, it is not normally used in the capture or intermediate steps of manufacturing processes. It is sometimes used as a final polishing step where a target protein is being separated from its aggregates or other significantly different molecular weight impurities. Another application area is the desalting of the purified target protein.


Tosoh Bioscience offers a number of Toyopearl HW-type products for size exclusion chromatography (Table I). These same SEC Toyopearl HW-type products are chemically modified with ion exchange groups, hydrophobic interaction or affinity ligands, to make the Toyopearl products shown in the later sections of this catalog.

Much of the information in this catalog section pertains to the available pore sizes, fractionation ranges, and particle sizes of the Toyopearl HW-type products, along with their physical and chemical properties. This data not only aids in the evaluation of the best resin for a SEC process step, but also lends insight into understanding the physical nature and the selection process of all Toyopearl resins.

FEATURES

- small particles available
- hydrophilic porous polymer structure
- narrow particle size distribution
- good mechanical stability
- chemically stable (pH 2 – 14)
- identical resin structure to TSKgel HPLC resins

FIGURE 1
RESIN CHEMISTRY OF TOYOPEARL SEC RESINS
(HYDROXYLATED ACRYLIC)



RESIN CHEMISTRY

Toyopearl size exclusion resins are highly hydroxylated polymethacrylic polymer beads (Figure 1). Their surface hydroxyl groups render them very hydrophilic and useful for protein separations. Toyopearl products including the functionalized materials seen in later catalog sections, have the least non-specific binding of any chromatographic resin. This is of particular note for separations such as blood factors where backbone interactions with the feedstock may result in decreased recovery of the targets. Their semi-rigid polymeric nature also gives them better pressure-flow characteristics than softer materials such as agarose.

PORE SIZE

Commercial Toyopearl HW-type size exclusion materials are available in 5 pore sizes covering 5 different fractionation ranges. The choice of Toyopearl HW products depends on the molecular weight of the feedstock components. Tables I and II show this information for proteins, dextrans and PEG polymers.

BENEFITS

- high resolution
- minimal non-specific adsorption effects
- high performance SEC – more efficient separations
- better pressure-flow characteristics
- excellent flow characteristics in large industrial size columns (up to 7 bar)
- constant packing volume over a wide range of salt concentrations
- compatible with organic solvents, can be cleaned in place (CIP) with acid or base
- stable polymer may be run at elevated temperature (4° - 60°C)
- autoclavable at 121 °C
- direct scale-up from TSKgel HPLC columns

SIZE EXCLUSION CHROMATOGRAPHY




TABLE I

PROPERTIES AND MOLECULAR WEIGHT SEPARATION RANGES FOR TOYOPEARL HW-TYPE RESINS (HW = HYDROPHILIC, WATER-COMPATIBLE POLYMERIC BASE RESINS)

Toyopearl resin	Particle size (µm)	Pore size (Å)	Molecular weight of sample (Da)		
			Polyethylene glycols and oxides	Dextrans	Globular proteins
HW-40S HW-40F HW-40C	20 - 40 30 - 60 50 - 100	50	100 - 3,000	100 - 7,000	100-10,000
HW-50S HW-50F	20 - 40 30 - 60	125	100 - 18,000	500 - 20,000	500 - 80,000
HW-55S HW-55F	20 - 40 30 - 60	500	100 - 150,000	1,000 - 200,000	1,000 - 700,000
HW-65S HW-65F HW-65C	20 - 40 30 - 60	1000	500 - 1,000,000	10,000 - 1,000,000	40,000 - 5,000,000
HW-75F	30 - 60	> 1000	4,000 - 5,000,000	100,000 - 10,000,000	500,000 - 50,000,000

TABLE II



a) GLOBULAR PROTEINS

* HW-55 is base material for some IEC and HIC products

** HW-65 is base material for most IEC, HIC and AFC products

CALIBRATION MOLECULES

B) DEXTRANS

C) POLYETHYLENE GLYCOLS



SIZE EXCLUSION CHROMATOGRAPHY

The Toyopearl HW-type resin range spans peptide and protein molecular weights between 100 - 50,000,000 Daltons. Each Toyopearl HW-type resin displays a typical calibration curve and exclusion limit for globular proteins (Figure 2).

Toyopearl HW-65 and HW-55 resins are the base beads for many Toyopearl products listed in the later sections of this catalog. A product having "-650" in its name uses the Toyopearl HW-65 resin for the base bead. Products having "-550" are derived from the Toyopearl HW-55 bead. The impact of pore size and its relationship to dynamic binding capacity and resin selectivity for a particular protein application is discussed in subsequent catalog sections.

See the hydrophobic interaction chromatography section of this catalog for Toyopearl PPG-600, Toyopearl Phenyl-600 and Toyopearl Butyl-600 resins which are functionalized on the Toyopearl HW-60 bead (not available for use in SEC).

PARTICLE SIZE

Resolution increases with decreasing particle size (Figure 3). Resin particle size is proportional to HETP and inversely proportional to the column efficiency and resolution of two peaks.

Most Toyopearl HW-type resins are available in three particle size ranges:

S-grade = 20 - 40 μm (Superfine)

F-grade = 30 - 60 μm (Fine)


C-grade = 50 - 100 μm (Coarse)

When the highest resolution is needed, the smaller S and F grade beads are preferred for process SEC. For desalting, where the resin is used in a filtration mode to remove the target from a buffer, the C grade is primarily employed because of its better flow dynamics at lower operating pressures.

Some Toyopearl HW-type products are also manufactured in "M-grade" (40-90 μm) and "EC-grade" (100-300 μm) to be used as the base beads for functionalized materials. These special grades are not commercially offered as SEC products. Please note that for the functionalized base beads in later catalog sections a "C-grade" is specified as a 50-150 μm bead and is not the SEC "C-grade" range of 50-100 μm .

FIGURE 2


CALIBRATION CURVES FOR GLOBULAR PROTEINS ON TOYOPEARL



Column: 22 mm ID x 30 cm L; Sample: Protein standards
Mobile phase: 0.06 mol/L phosphate buffer (pH 7.0) in 0.06 mol/L KCl;
Legend: V_e = elution volume, V_c = column volume;
Detection: UV @ 280 nm

FIGURE 3

COMPARISON OF RESOLUTION ON DIFFERENT PARTICLE SIZES OF TOYOPEARL HW-55



Column: Toyopearl HW-55, 26 mm ID x 70 cm L
Sample: 1. Thyroglobulin (0.3%), 2. γ -Globulin (0.3%),
3. β -Lactoglobulin (0.3%), 4. Cytochrome C (0.1%)
Mobile phase: 33.3 mmol/L phosphate buffer (pH 7.0), 0.2 mol/L NaCl
Flow rate: 106 mL/h (20 cm/h); Inj. vol.: 1 mL;
Temperature: 25°C; Detection: UV @ 280 nm

SIZE EXCLUSION CHROMATOGRAPHY




MOBILE PHASE

Mobile phase components, such as salts, can affect SEC separations. The presence or absence of sodium chloride influences the elution volume of proteins. This is demonstrated in Figure 4, in which a mixture of various proteins was separated on a column packed with Toyopearl HW-55F. Salt concentrations can change the hydrodynamic radius of proteins and either increase or decrease their molecular size as a function of salt strength. Ideally, in SEC sample components do not interact with the packing material. In practice it is often necessary to select a salt concentration which minimizes secondary interactions of the sample components with the resin. However, there are instances where secondary interactions, particularly hydrophobic interactions at higher salt concentrations, can be exploited.

It is important to note that relatively minor changes in protein structure may affect protein solubility and encourage secondary hydrophobic interactions causing similarly sized proteins or analogs to elute at different times. In those cases it may be necessary to modify the mobile phase composition to regain a separation based on molecular size alone.

FIGURE 4

COMPARISON OF THE ELUTION VOLUMES OF PROTEINS IN PRESENCE AND ABSENCE OF NaCl



Column: Toyopearl HW-55F, 22 mm ID x 50 cm L; Elution: 25 mmol/L Tris-HCl with or without 0.5 mol/L NaCl, (pH 7.5); Flow rate: 16 cm/h; Temperature: 5 - 10 °C; Detection: UV @ 280 nm, 420 nm for heme proteins, 200 nm for proteins without aromatic amino acid

PROPERTIES OF TOYOPEARL SEC RESINS IN AQUEOUS ELUENTS

- **HIGH MECHANICAL STABILITY**
Toyopearl resins can be operated at pressures up to 3 bar without deformation.
- **MINIMUM CHANGE IN GEL BED VOLUME**
Changes in the column bed volume under operational salt conditions are negligible. Toyopearl does not shrink or swell even in high concentrations of strong denaturing agents such as urea or guanidine hydrochloride.
- **CHEMICAL STABILITY**
Toyopearl is stable from pH 2-13, and tolerant to pH 0-14 for short periods. Biomolecules which are only soluble at extreme pH values can be readily separated.
- **SHARP CHROMATOGRAPHIC PEAKS**
Toyopearl's narrow particle size distribution (min. 80% – within declared limits) results in better peak shapes and higher elution target concentrations than other SEC materials.
- **TEMPERATURE STABILITY**
Toyopearl is thermally stable and does not degrade or denature even in boiling water. Toyopearl resins can be sterilized by autoclaving at 121°C.
- **MICROORGANISM RESISTANCE**
Toyopearl is an organosynthetic material, and is resistant to degradation by microorganisms.
- **SUITABILITY FOR ENZYME IMMOBILIZATION**
Toyopearl resins contain numerous hydroxyl groups on the external and internal bead surfaces. These, in combination with the chemical stability of the polymer, make the resin well suited for the covalent bonding of enzymes or other ligands. (Please see the AFC section for more information.)



SIZE EXCLUSION CHROMATOGRAPHY

PROPERTIES IN ORGANIC ELUENTS

Toyopearl resins can be used in organic solvents or mixtures of organic solvents and water. Bed volumes may swell or shrink relative to water depending on the solvent as shown in Tables III and IV. DMSO can be used for SEC of oligosaccharides and polyethylene glycols. The compatibility of DMF with Toyopearl also permits SEC separation of hydrophobic substances such as polystyrenes.

All of the physical and chemical properties listed above for the Toyopearl HW-type series of SEC products makes them an excellent choice for use as the base beads for the ion exchange, hydrophobic interaction, and affinity chromatographic resins discussed in the later sections of this catalog.

OTHER APPLICATIONS

The Toyopearl HW-type resins are commonly used in size exclusion chromatography and desalting applications. Some other important uses of these materials are:

- Removal of surfactants such as Triton® X-100 from biological solutions by an adsorption mechanism
- Use in hydrophobic interaction chromatography (HIC) for the separation of very hydrophobic molecules
- Use in HIC separations as a guard column for hydrophobic impurities
- Possible use as a stationary phase for either normal or reversed phase separations depending on solvent system selected

➤ **TABLE III**

SWELLING PROPERTIES IN VARIOUS SOLVENTS

Toyopearl	HW-40	HW-50	HW-55	HW-65	HW-75
Water	100	100	100	100	100
0.2 mol/L KCl	100	100	100	100	100
MeOH	100	100	100	100	105
EtOH	100	100	100	100	110
DMF	110	110	105	105	120
Acetone	80	80	85	90	110
Toluene	65	70	70	75	90

➤ **TABLE IV**

ADDITIONAL SWELLING DATA FOR TOYOPEARL HW-40

Toyopearl	DMSO	Ethyl Acetate	Benzene	CHCl ₃	CHCl ₃ /MeOH (1:1)
HW-40	140	80	70	105	120

SIZE EXCLUSION CHROMATOGRAPHY



ORDERING INFORMATION


TOYOPEARL LABPAK

PART #	PRODUCT DESCRIPTION	CONTAINER SIZE (mL)	PARTICLE SIZE (µm)
19821	SECPAK LMW (HW-40F, HW-50F, HW-55F)	3 x 150	30 - 60
19819	SECPAK HMW (HW-55F, HW-65F, HW-75F)	3 x 150	30 - 60
19820	SECPAK HP (HW-40S, HW-50S, HW-55S, HW-65S)	4 x 150	20 - 40

TOYOPEARL SEC RESINS:

Conditions: Exclusion limits are +/- 30% and are determined using PEG, PEO, or dextran standards, as appropriate.

PART #	PRODUCT DESCRIPTION	CONTAINER SIZE (mL)	PARTICLE SIZE (µm)	EXCLUSION LIMIT (Da)
19809	Toyopearl HW-40S	150	20 - 40	3 x 10 ³
07451		250		
14681		1,000		
07967		5,000		
19808	Toyopearl HW-40F	150	30 - 60	3 x 10 ³
07448		500		
14682		1,000		
07968		5,000		
19807	Toyopearl HW-40C	150	50 - 100	3 x 10 ³
07449		500		
14683		1,000		
07969		5,000		
19811	Toyopearl HW-50S	150	20 - 40	1.8 x 10 ⁴
07455		250		
14684		1,000		
08059		5,000		
19810	Toyopearl HW-50F	150	30 - 60	1.8 x 10 ⁴
07453		500		
14685		1,000		
08060		5,000		
18368		50,000		
19813	Toyopearl HW-55S	150	20 - 40	1.5 x 10 ⁵
07459		250		
14686		1,000		
08062		5,000		
19812	Toyopearl HW-55F	150	30 - 60	1.5 x 10 ⁵
07457		500		
14687		1,000		
08063		5,000		
19815	Toyopearl HW-65S	150	20 - 40	1 x 10 ⁶
07467		250		
14688		1,000		
08068		5,000		
19814	Toyopearl HW-65F	150	30 - 60	1 x 10 ⁵
07465		500		
14689		1,000		
08069		5,000		
21481	Toyopearl HW-65C	150	50 - 100	1 x 10 ⁵
07466		500		
14690		1,000		
08070		5,000		
21482		50,000		
19816	Toyopearl HW-75F	150	30 - 60	8.25 x 10 ⁵
07469		500		
14691		1,000		
08072		5,000		



IEC ION EXCHANGE CHROMATOGRAPHY

IEC PRODUCTS

➤ ANION EXCHANGE

- Toyopearl SuperQ-650
- Toyopearl QAE-550
- Toyopearl Q-600C AR
- Toyopearl DEAE-650
- Toyopearl GigaCap Q-650
- TSKgel SuperQ-5PW
- TSKgel DEAE-5PW


➤ CATION EXCHANGE

- Toyopearl MegaCap II SP-550EC
- Toyopearl SP-650
- Toyopearl SP-550
- Toyopearl CM-650
- Toyopearl GigaCap S-650
- Toyopearl GigaCap CM-650
- TSKgel SP-3PW
- TSKgel SP-5PW

≡ TOSOH FACT

The Tosoh logo symbolizes the corporate philosophy of Tosoh's vision of the ideal.

The curved lines represent the realization of happiness, reflecting Tosoh's management philosophy of putting people first. The square in the center expresses the advanced nature of Tosoh's technology and also represents the outstanding quality of Tosoh's products. The right-angle cut at the top portrays an image of contributing to society, Tosoh's stance towards the outside world. The red corporate color symbolizes the Tosoh spirit, which guides the ceaseless efforts to realize the ideal.





ION EXCHANGE CHROMATOGRAPHY

TOYOPEARL ION EXCHANGE CHROMATOGRAPHY RESINS

Ion Exchange Chromatography (IEC) is the most common liquid chromatographic method used in manufacturing therapeutic proteins. Due to the high dynamic binding capacities of ion exchange resins relative to those of the other chromatographic modes (Table I), it is the chromatographic technique selected by many developers for the capture or concentration step. Tosoh Bioscience offers a broad range of products for ion exchange applications.

HOW DOES IEC WORK?

IEC is based on the binding of proteins to positively or negatively charged groups which are immobilized on a stationary phase and which are in equilibrium with free counter ions in the mobile phase. In the process of adsorption, the mobile phase counter ions are exchanged by the protein solute. The binding of proteins to the ion exchange matrix predominantly occurs via charged amino acid residues located at the surface of the protein molecule.

The development of optimum chromatographic system conditions requires knowledge of both the protein's pI and the pKa of the ion exchange media. A binding buffer pH is selected between the pI of the target and the ion exchanger's pKa (Figure 1). This ensures that the protein is in the opposite charge state compared to the ion exchange media. When possible, the pH is also optimized to effect the highest solubility of the target protein. Higher protein solubilities make more efficient use of the overall ion exchange capacity of the resin. A salt is selected as the source of counter ions in the mobile phase and elution occurs as the salt strength is increased to a higher concentration than the target's binding salt conditions.

ION EXCHANGE GROUPS AVAILABLE

Toyopearl and TSKgel 5PW-type IEC resins are available with six different ion exchange groups as shown in Table II:

- 3 for anion exchange – Q, QAE, DEAE
- 3 for cation exchange – S, SP, CM

FEATURES

- porous, hydrophilic polymer based resin
- chemical stability
- column bed stability
- mechanical stability
- continuous selectivity

FIGURE 1

PK_a VALUES FOR ION EXCHANGE GROUPS




TABLE I

Separation Mode	Binding capacity for standard proteins (mg/mL)	Binding capacity in production processes (mg/mL)
Ion Exchange	100-200	50-100
Hydrophobic Interaction	40-60	10-30
Affinity (group specific ligands)	40-100	20-60
Reversed Phase (polymeric media)	60-100	30-50

PORE SIZES

Tosoh Bioscience offers a range of pore sizes for IEC resins originating from our size exclusion chromatography base resins. Four different mean pore diameters are used for the current ion exchange resins: 1000 Å, 750 Å, 500 Å, and 250 Å (Table III). Depending on the kind of ligand attachment, the effective pore size of the resulting IEC resin is smaller than the pore size of the base bead. When network ligand technology is applied the accessible pore size is varying with pH and salt concentration, therefore all pore sizes mentioned here are those of the respective Toyopearl HW or TSKgel base resin.

BENEFITS

- suitable for laboratory scale and process chromatography
- autoclavable at 121°C
- temperature range 4° - 60°C
- pH range 2-13, can be regenerated with acid or base
- compatible with organic solvents
- constant packing volume over a wide range of salt concentrations
- excellent flow characteristics in large industrial columns (up to 7 bar)
- easy scale-up from TSKgel IEC columns
- high yields of biologically active proteins

ION EXCHANGE CHROMATOGRAPHY


TABLE II

ANION EXCHANGE RESIN	Base bead	Pore size (Å)	Particle size (µm)	IEC (meq/mL)	DBC (BSA) (mg/mL-gel)	Recovery (%)
Toyopearl GigaCap Q-650M	HW-65	1,000	50-100	0.17	172	97
Toyopearl SuperQ-650M	HW-65	1,000	40-90	0.24	145	98
Toyopearl DEAE-650M	HW-65	1,000	40-90	0.11	25	97
Toyopearl Q-600C AR	HW-60	750	50-150	0.18	100	98
Toyopearl QAE-550C	HW-55	500	50-150	0.36	30	96
CATION EXCHANGE RESIN	Base bead	Pore size (Å)	Particle size (µm)	IEC (meq/mL)	DBC (hIgG) (mg/mL-gel)	Recovery (%)
Toyopearl GigaCap S-650M	HW-65	1,000	50-100	0.16	145	98
Toyopearl GigaCap CM-650M	HW-65	1,000	50-100	0.24	100	98
Toyopearl SP-650C	HW-65	1,000	50-150	0.12	12	98
Toyopearl CM-650C	HW-65	1,000	50-150	0.09	35 (lysozyme)	
Toyopearl SP-550C	HW-55	500	50-150	0.13	14	98

TABLE III

Base bead	Toyopearl HW-65 TSKgel G5000PW	Toyopearl HW-60	Toyopearl HW-55	TSKgel G3000PW
Pore Diameter	1000 Å	750 Å	500 Å	250 Å
	Toyopearl GigaCap S-650 Toyopearl GigaCap CM-650 Toyopearl GigaCap Q-650 Toyopearl SuperQ-650 Toyopearl DEAE-650 Toyopearl SP-650 Toyopearl CM-650 TSKgel SP-5PW TSKgel SuperQ-5PW TSKgel DEAE-5PW	Toyopearl Q-600C AR	Toyopearl SP-550 Toyopearl MegaCap II SP-550 Toyopearl QAE-550	TSKgel SP3-PW

HIGHER ACCESSIBLE SURFACE AREA EQUALS MORE CAPACITY


A bead with a small pore size has theoretically more surface area than the same size bead with a larger pore. Figure 2 shows insulin binding capacity on six different pore size beads. As the pore size increases to the point where the insulin has maximal access to the internal surface area the insulin capacity increases.

However, there is a point of diminishing return. Because the absolute surface area decreases as the pores become larger, the insulin capacity decreases accordingly. This effect is similarly demonstrated by the lysozyme static binding capacity range specifications of Toyopearl SP-550C and Toyopearl SP-650C which are 80-120 g/L and 35-55 g/L, respectively (data not shown).

For more information on pore size and particle size combinations not contained in our commercial products, please see the Custom Resin section of this catalog.

FIGURE 2

OPTIMIZATION OF INSULIN BINDING CAPACITY AS A FUNCTION OF PORE SIZE OF EXPERIMENTAL TSKgel SP-TYPE RESINS





ION EXCHANGE CHROMATOGRAPHY

MULTIPLE PARTICLE SIZES SIMPLIFY SCALING UP OR DOWN

Because Toyopearl HW-65 and TSKgel G5000PW products have similar backbone chemistry and selectivity, scaling up or scaling down for a selected ion exchange method is simple. Practically speaking, the chromatographic conditions that work for one particle size will work for all particle sizes with a given ligand functionality. The elution order of the feedstock components will remain the same with increasing resolution as the particle size gets smaller (Figure 3). Figure 4 lists the complete range of ion exchange products, particle sizes and suggests how they are typically placed into a manufacturing scheme. The availability of smaller bead sizes for greater resolution while maintaining the same selectivity is particularly useful in the areas of oligonucleotide and peptide purification.

RESIN PHYSICAL PROPERTY SELECTION

(for resins available in different pore sizes with the same ligand and ligand attachment chemistry)

For bind/elute chromatography:

- Select the smallest pore size resin appropriate for the size of the target molecule.
- Select a larger particle size for a capture step, a smaller one for intermediate or polishing steps

For flow through chromatography:

- If the target molecule's size is larger than most components of the feedstream, select a pore size which will tend to exclude it (known as kinetic exclusion, this technique can also be used under binding conditions as the excluded molecule only sees 1% of the resin surface area and the capacity/recovery loss is minimal)

For large molecule impurity clearance:

- Select a pore size which includes the target molecule, but excludes the impurity. (see the calibration curves of the Toyopearl base beads in the SEC section of the catalog as an aid).

OLIGONUCLEOTIDE PURIFICATIONS


Table 4 shows the different particle sizes available with the anion exchange SuperQ functionality, which is typically used for oligonucleotide purifications. The relative binding capacities and predicted resolution of the five particle sizes are depicted by a series of "+" characters. The more "+" characters listed in the table the better one resin is relative to another for that parameter. If a process is developed using one of the five resins and more resolution is needed, select an appropriate smaller particle size product. Similarly if more capacity is needed, and resolution is not a critical issue, a

larger particle size resin can be selected. The larger particle Toyopearl resins are less crosslinked than the corresponding TSKgel 5PW type resin products and have more active sites for ligand attachment. Thus they have higher capacities than the TSKgel 5PW-type resins. In some cases, Toyopearl GigaCap Q-650M (also shown in Table IV) and its very high capacity can be used, although its selectivity is somewhat different than the other quaternary Q products.

TSKgel SuperQ-5PW products typically have 2-4 times the binding capacity of other small particle anion exchange resins available on the market. This has significant bearing in the area of difficult to resolve "n-1" DNA and RNAi purifications as loading amounts are increased. Under higher loading conditions, the TSKgel SuperQ resins maintain their resolution much better than smaller particle, lower capacity resins. The smaller particle products may start out with a slight separation advantage under low oligonucleotide loading conditions, but this vanishes as the feedstock load is increased.

Figure 5 shows a comparison of one smaller particle size, competitive product, which starts out having better resolution than TSKgel SuperQ-5PW (20) resin at 1 mg oligonucleotide/mL of resin. At 20 mg oligonucleotide/mL of resin, however, the resolution of peaks on the competitive product deteriorates significantly. The TSKgel SuperQ-5PW (20) retains excellent resolution even at this higher oligonucleotide level.

FIGURE 3 COMPARISON OF TSKgel DEAE-5PW AND TOYOPEARL DEAE-650S



Column: 55 mm ID x 20 cm L; Sample: calf liver acetone powder, 94 mg in 4.7 mL in 0.02 mol/L Tris-HCl (pH 8.0)
 Mobile phase: 100 min linear gradient from 0 mol/L to 0.25 mol/L NaCl followed by a 20 min linear gradient from 0.25 mol/L to 0.5 mol/L NaCl in 0.02 mol/L Tris-HCl (pH = 8.0)
 Flow rate: 50 cm/h; Detection: UV @ 280 nm

ION EXCHANGE CHROMATOGRAPHY



FIGURE 4

Process step	Bead size (µm)	Process media	
		ANION	CATION
Capture	200 µm		Toyopearl MegaCap II SP-550EC
	100 µm	Toyopearl SuperQ-650C Toyopearl DEAE-650C Toyopearl QAE-550C	Toyopearl SP-650C Toyopearl SP-550C Toyopearl CM-650C
	75 µm	Toyopearl GigaCap Q-650M	Toyopearl GigaCap S-650M
Intermediate Purification	65 µm	Toyopearl SuperQ-650M Toyopearl DEAE-650M	Toyopearl SP-650M Toyopearl CM-650M
	35 µm	Toyopearl SuperQ-650S Toyopearl DEAE-650S	Toyopearl SP-650S Toyopearl CM-650S
	30 µm	TSKgel SuperQ-5PW (30) TSKgel DEAE-5PW (30)	TSKgel SP-5PW (30) TSKgel SP-3PW (30)
Polishing	20 µm	TSKgel SuperQ-5PW (20) TSKgel DEAE-5PW (20)	TSKgel SP-5PW (20)
	10 µm	TSKgel SuperQ-5PW 7.5 mm ID x 7.5 cm L TSKgel DEAE-5PW 7.5 mm ID x 7.5 cm L	TSKgel SP-5PW 7.5 mm ID x 7.5 cm L TSKgel CM-5PW 7.5 mm ID x 7.5 cm L

Same selectivity HPLC columns are available for most process media

TABLE IV

Oligonucleotide Purification Products:

	Bead size (mean µm)	Binding capacity	Resolution	Bead type	Ligand attachment
TSKgel SuperQ-5PW (20)	20	++	+++++	methacrylic	Type A
TSKgel SuperQ-5PW (30)	30	++	++++	methacrylic	Type A
Toyopearl SuperQ-650S	35	++++	+++	methacrylic	Type A
Toyopearl SuperQ-650M	65	++++	++	methacrylic	Type A
Toyopearl SuperQ-650C	100	++++	+	methacrylic	Type A
Toyopearl GigaCap Q-650M	75	+++++	++	methacrylic	Type B

Peptide Purification Products:

	Bead size (mean µm)	Binding capacity	Resolution	Bead type	Ligand attachment
TSKgel SP-5PW (20)	20	++	+++++	methacrylic	Traditional
TSKgel SP-5PW (30)	30	++	++++	methacrylic	Traditional
TSKgel SP-3PW (30)	30	+++	++++	methacrylic	Traditional
Toyopearl SP-650S	35	++++	+++	methacrylic	Traditional
Toyopearl SP-650M	65	++++	++	methacrylic	Traditional
Toyopearl SP-650C	100	++++	+	methacrylic	Traditional
Toyopearl GigaCap S-650M	75	+++++	++	methacrylic	Type B



ION EXCHANGE CHROMATOGRAPHY

PEPTIDE PURIFICATIONS

Cation exchange chromatography is commonly used for peptide purification. Table 4 shows the same particle size availability of Toyopearl and TSKgel PW strong cation exchange resins. Based on the needs for capacity and resolution, an appropriate S or SP resin is selected for a particular peptide application. The new TSKgel SP-3PW (30) is based on a 250 Å pore resin. It was developed to provide high binding capacities for peptides and small proteins. It also has a different selectivity than TSKgel SP-5PW (30).

LIGAND ATTACHMENT TECHNOLOGY

Tosoh Bioscience applies three different ligand attachment chemistries reflecting three generations of ligand attachment technology. The “traditional” method, or first generation of attaching the ion exchange ligand, is directly to the resin surface through a proprietary spacer arm. Toyopearl and TSKgel PW type ion exchange resins using this traditional bead functionalization method are:


- Toyopearl SP-650 and SP-550
- TSKgel SP-3PW and SP-5PW
- Toyopearl CM-650
- Toyopearl Q-600C AR
- Toyopearl Q-550
- Toyopearl DEAE-650
- TSKgel DEAE-5PW
- Toyopearl MegaCap II SP-550EC

A second generation attachment chemistry (Type A noted in Table IV) for increasing protein binding within the accessible surface area, is to add a carbon spacer network between the bead surface and the ligand. It is also possible to attach ligand groups along the length of the spacer network thus improving capacity. There are two resins which incorporate this type A ligand attachment chemistry:

- Toyopearl SuperQ-650
- TSKgel SuperQ-5PW

➤ FIGURE 5

TSKgel SuperQ-5PW (20) MAINTAINS RESOLUTION AT HIGH OLIGONUCLEOTIDE LOAD



Column: 0.66 cm x 15 cm L (5.1 mL) (resin as noted in figure); Flow rate: 1.43 mL/min (250 cm/h)

Buffer A: 20 mmol/L Tris-HCl + 10 mmol/L EDTA (pH= 9.0); Buffer B: 20 mmol/L Tris-HCl + 10 mmol/L EDTA + 1.0 mol/L NaCl (pH= 9.0)

Sample loaded: DNA based oligonucleotides were loaded as followed: 1 mg/column panels A & B, 20 mg/column panels C & D

Separation conditions: Column was washed with 5 CV 100% buffer A followed by 11 mL injection. Column was then washed with 3 CV 100% buffer A followed by 6 CV of linear gradient 35-53 buffer B. Finally, column was washed with 5 CV 100% buffer B.

Detection: UV @ 254 nm; Fractions: 0.5 mL fractions were taken from peaks of interest and analyzed on a TSKgel DNA-NPR column

ION EXCHANGE CHROMATOGRAPHY



A third generation ligand attachment chemistry (Type B noted in Table IV) improves the accessible location of the ligand groups. The result of this modification is significantly increased capacity and improved mass transfer. Improved mass transfer also reduces the target molecule elution volume. All Toyopearl GigaCap resins use this Type B ligand attachment chemistry:

- Toyopearl GigaCap S-650
- Toyopearl GigaCap CM-650
- Toyopearl GigaCap Q-650


Figures 6, 7, and 8 show the breakthrough curves for the three Toyopearl GigaCap resins. They are compared where possible with the most current equivalent competitive resin. Each trace shows the dynamic binding capacity of the resin up to 10% breakthrough plus the elution profile for the target molecule. Please note the significant reduction in elution pool volumes of the Toyopearl GigaCap resins when compared to other products. The concentration of the eluted peak is proportionally increased as well. It is possible to achieve reductions in elution pool volumes of over 75%. This can reduce the cost of further downstream process steps.

IMPACT OF PORE SIZE AND LIGAND ATTACHMENT ON DBC

Table 5 contains DBC data for five of our Toyopearl resins using three different sized proteins. There are three different pore sizes and three different ligand attachment methods represented. Toyopearl GigaCap Q-650M has the highest capacity for all combinations of pore size and attachment chemistries. Please note the decrease in capacity for the larger proteins on the Toyopearl SuperQ-650M resin indicating that the accessible pore volume has diminished by the ligand attachment chemistry used.

➤ **FIGURE 6**

TOYOPEARL GIGACAP S-650M VS. COMPETITOR S ELUTION POOL VOLUME COMPARISON



Column size: 6 mm ID x 40 mm bed; Sample: polyclonal human IgG (1 mg/mL); Loading Buffer: 0.1 mol/L acetate buffer (pH= 4.7) Elution Buffer: 0.1 mol/L acetate buffer (pH= 4.7) + 1.0 mol/L NaCl Linear velocity: 212 cm/h; Detection: UV @ 280 nm

PEGylated PROTEINS


Ion exchange resins are frequently used for the purification of pegylated proteins. Figure 9 shows the breakthrough curves of five Toyopearl cation exchange resins for mono-pegylated lysozyme.

RESIN PRESSURE FLOW PROPERTIES

All Toyopearl resins are designed to withstand pressures up to 3 bar. The newer Toyopearl GigaCap resins have a particle size of 50-100 microns which is slightly larger than our normal M-grade 40-90 micron beads. This particle size difference generates a lower back pressure (Figure 10) than our more traditional M-grade ion exchange products. The TSKgel 5PW type resins can be operated at pressures up to 20 bar.

➤ **FIGURE 7**


TOYOPEARL GIGACAP Q-650M VS. COMPETITOR Q ELUTION POOL VOLUME COMPARISON



Column size: 6 mm ID x 40 mm L; Sample: polyclonal human IgG (1 mg/mL); Loading buffer: 15 mmol/L Tris-HCl (pH= 8.7); Elution buffer: 15 mmol/L Tris-HCl (pH= 8.7) + 1.0 mol/L NaCl; Linear velocity: 212 cm/h; Detection: UV @ 280 nm

➤ **FIGURE 8**

TOYOPEARL GIGACAP CM-650M ELUTION POOL VOLUME



Column size: 6 mm ID x 40 mm L; Sample: polyclonal human IgG (1 mg/mL); Loading buffer: 50 mmol/L sodium acetate buffer (pH= 4.7); Elution buffer: 15 mmol/L Tris-HCl (pH= 8.7) + 1.0 mol/L NaCl Linear velocity: 212 cm/h; Detection: UV @ 280 nm



ION EXCHANGE CHROMATOGRAPHY


If recommended packing procedures are followed, Toyopearl and TSKgel IEC resins maintain stable bed volumes during the pH and ionic strength changes that occur during normal ion exchange chromatography (Consult our Toyopearl Instruction Manual for the recommended packing conditions). Multi-cycle gradient operation and re-equilibration are accomplished without volume changes in the Toyopearl column bed. The mechanical stability of Toyopearl resins allows the use of longer column beds with more efficiency or higher operational flow rates.

ALKALINE STABILITY

Tosoh has focused on improving the alkaline stability of its newer ion exchange resins. Higher capacity resins can bind not only more of the target molecule, but the impurities and isoforms as well. In some cases more rigorous cleaning agents like 0.5 mol/L NaOH and even 1.0 mol/L NaOH are needed to insure proper resin regeneration. Naturally, the resins need to tolerate these more stringent conditions. As seen in Table VI all three of the new Toyopearl GigaCap series ion exchange resins have excellent alkaline stability.

FIGURE 9


BREAKTHROUGH CURVES OF MONO-PEGYLATED LYSOZYME USING TOYOPEARL CATION EXCHANGE RESINS



Dynamic binding capacities were determined at 10% breakthrough. Column size: 6 mm ID x 40 mm L; Sample: mono-PEGylated lysozyme; Loading buffer: 20 mmol/L phosphate buffer (pH= 7.0) Elution buffer: 20 mmol/L phosphate buffer (pH= 7.0) + 0.5mol/L NaCl. Linear velocity: 212 cm/h; Detection: UV @ 280 nm PEG MW= 5kDa

FIGURE 10

PRESSURE-FLOW CURVE COMPARISON



Column size: 22 mm ID x 20 cm L; Mobile phase: distilled water. Temperature: 25 °C

TABLE V

DYNAMIC BINDING CAPACITY VARIES WITH PROTEIN SIZE

Resin	Pore size (Å)	BSA 66 kDa	Binding capacity (mg/mL gel)	
			Human IgG 160 kDa	Thyroglobulin 660 kDa
Toyopearl GigaCap Q-650M	1.000	173	108	71
Toyopearl SuperQ-650M	1.000	145	13	3
Toyopearl Q-600C AR	750	108	90	26
Toyopearl QAE-550C	500	29	32	6
Toyopearl DEAE-650M	1.000	25	31	3

Column size: 6 mm ID x 4 cm L; Sample concentration: 1 mg/mL; Loading buffers: BSA 0.05 mol/L Tris-HCl (pH = 8.5); Human IgG 0.05 mol/L Tris-HCl (pH = 8.7); Thyroglobulin 0.05 mol/L Tris-HCl (pH = 8.7) + 0.15 mol/L NaCl

Elution buffers: loading buffer + 1.0 mol/L NaCl; Flow rate: 212 cm/h; Detection: UV @ 280 nm

ION EXCHANGE CHROMATOGRAPHY



TABLE VI

TOYOPEARL GIGACAP RESIN BASE STABILITY

Resin solution	Storage Molecule	Test	Capacity	Starting capacity	Week 1	Week 2	Week 3
Toyopearl GigaCap S-650M	1.0 mol/L NaOH	h-IgG	Dynamic	143 (mg/mL)	144	140	135
Toyopearl GigaCap CM-650M	0.5 mol/L NaOH	h-IgG	Dynamic	99 (mg/mL)	88	90	91
Toyopearl GigaCap Q-650M	0.5 mol/L NaOH	BSA	Static	166 (mg/mL)	NA	153*	136


* 12 days

FIGURE 11

TOYOPEARL Q-600C AR

TOYOPEARL Q-600C AR RESIN DBC AS A FUNCTION OF SODIUM HYDROXIDE EXPOSURE

For this catalog we are announcing a new, high capacity anion exchange resin, Toyopearl Q-600C AR (using first generation ligand attachment chemistry). This new resin has a slightly smaller pore size than Toyopearl GigaCap Q-650M and has a typical BSA binding capacity of 100 mg/mL. As shown in Figure 11, after 100 days of exposure to 1 mol/L NaOH, its DBC remains unchanged. Figure 12 shows the preservation of selectivity after extensive exposure to caustic.



TOYOSCREEN PREPACKED COLUMNS FOR PROCESS DEVELOPMENT


ToyoScreen columns packed with the full range of our Toyopearl IEC products are available in 1 mL and 5 mL resin volumes. The ToyoScreen columns provide a convenient way to perform early resin screening for both target retention and recovery. Multiple columns can be connected in series for additional capacity or resolution. Please see the ordering information at the end of this section or contact us for more information on these products.

FIGURE 12

LABPAK

STABILITY OF TOYOPEARL Q-600C AR RESIN AFTER EXPOSURE TO 1 MOL/L NaOH

For scientists wishing to develop a better physical understanding of the packing properties of Toyopearl and TSKgel ion exchange resins, we offer Toyopearl LabPaks containing small volumes of the bulk resins. Please see the ordering information at the end of this section or contact us for more information on these products.



SUMMARY

Since 2007, major additions have been made to the Toyopearl ion exchange resins offered by Tosoh Bioscience. The Toyopearl GigaGap family of resins and the new Toyopearl Q-600C AR media represent significant improvements in terms of dynamic binding capacity, elution kinetics, and alkaline stability. When these products are used, process developers have more options available to design and optimize their process to improve productivity and lower operating costs.

Column: 6.0 mm ID x 4 cm L; Flow rate: 1.0 mL/min;
 Elution: Buffer A: 0.05 mol/L Tris-HCl buffer (pH= 8.5); Buffer B: 0.05 mol/L Tris-HCl buffer + 1.0 mol/L NaCl (pH= 8.5);
 Gradient: 60-min linear gradient from buffer A to buffer B;
 Detection: UV @ 280 nm



ION EXCHANGE CHROMATOGRAPHY

ORDERING INFORMATION

TOYOSCREEN PROCESS DEVELOPMENT COLUMNS FOR IEC:

PART #	PRODUCT DESCRIPTION	PACKAGE
21360	ToyoScreen DEAE-650M	1 mL x 6 each
21361	ToyoScreen DEAE-650M	5 mL x 6 each
21859	ToyoScreen GigaCap Q-650M	1 mL x 6 each
21860	ToyoScreen GigaCap Q-650M	5 mL x 6 each
21362	ToyoScreen SuperQ-650M	1 mL x 6 each
21363	ToyoScreen SuperQ-650M	5 mL x 6 each
21364	ToyoScreen QAE-550C	1 mL x 6 each
21365	ToyoScreen QAE-550C	5 mL x 6 each
21992	ToyoScreen Q-600C AR	1 mL x 6 each
21993	ToyoScreen Q-600C AR	5 mL x 6 each
21951	ToyoScreen GigaCap CM-650M	1 mL x 6 each
21952	ToyoScreen GigaCap CM-650M	5 mL x 6 each
21366	ToyoScreen CM-650M	1 mL x 6 each
21367	ToyoScreen CM-650M	5 mL x 6 each
21868	ToyoScreen GigaCap S-650M	1 mL x 6 each
21869	ToyoScreen GigaCap S-650M	5 mL x 6 each
21870	ToyoScreen MegaCap II SP-550EC	1 mL x 6 each
21871	ToyoScreen MegaCap II SP-550EC	5 mL x 6 each
21368	ToyoScreen SP-650M	1 mL x 6 each
21369	ToyoScreen SP-650M	5 mL x 6 each
21370	ToyoScreen SP-550C	1 mL x 6 each
21371	ToyoScreen SP-550C	5 mL x 6 each
21392	ToyoScreen IEC Anion Mix Pack, NEW (DEAE-650M, SuperQ-650M, QAE-550C, GigaCap Q-650M, Q-600C)	1 mL x 5 Grades x 1 each
21393	ToyoScreen IEC Anion Mix Pack (DEAE-650M, SuperQ-650M, QAE-550C, GigaCap Q-650M, Q-600C)	5 mL x 5 Grades x 1 each
21394	ToyoScreen IEC Cation Mix Pack NEW (CM-650M, SP-650M, SP-550C, GigaCap CM-650M, GigaCap S-650M)	1 mL x 5 Grades x 1 each
21395	ToyoScreen IEC Cation Mix Pack (CM-650M, SP-650M, SP-550C, GigaCap CM-650M, GigaCap S-650M)	5 mL x 5 Grades x 1 each
21396	ToyoScreen IEC Mix Pack, NEW (GigaCap Q-650M, GigaCap S-650M, GigaCap CM-650M, Q-600C, SuperQ-650M, SP-550C)	1 mL x 6 Grades x 1 each
21397	ToyoScreen IEC Mix Pack, NEW (GigaCap Q-650M, GigaCap S-650M, GigaCap CM-650M, Q-600C, SuperQ-650M, SP-550C)	5 mL x 6 Grades x 1 each

TOYOSCREEN COLUMN ACCESSORIES

PART #	PRODUCT DESCRIPTION
21400	ToyoScreen Column Holder



ION EXCHANGE CHROMATOGRAPHY



ORDERING INFORMATION

TSKgel LABPAK

PART #	PRODUCT DESCRIPTION	CONTAINER SIZE (mL)	PARTICLE SIZE (µm)
43380	IEXPAK PW (20) (SP-5PW, DEAE-5PW, SuperQ-5PW)	3 x 25 mL	15-25
43280	IEXPAK PW (30) (SP-5PW, DEAE-5PW, SuperQ-5PW)	3 x 25 mL	20-40

TOYOPEARL LABPAK

PART #	PRODUCT DESCRIPTION	CONTAINER SIZE (mL)	PARTICLE SIZE (µm)
19817	IEXPAK HP (CM-650S, SP-650S, DEAE-650S, SuperQ-650S)	4 x 25 mL	20-50
43210	AIEXPAK (GigaCap Q-650M, SuperQ-650M, Q-600C) NEW	3 x 100 mL	40-90 and 50-150
43220	CIEXPAK (GigaCap CM-650M GigaCap S-650M, SP-650M, SP-550C) NEW	3 x 100 mL	40-90 and 50-150

ANION EXCHANGE RESINS:

TOYOPEARL BULK MEDIA

PART #	PRODUCT DESCRIPTION	CONTAINER SIZE (mL)	PARTICLE SIZE (µm)	ION EXCHANGE CAPACITY (meq/mL RESIN)	TYPICAL CAPACITY (mg BSA/mL RESIN)
43271	Toyopearl QAE-550C	100	50-150	0.28-0.38	60-80
14026		250			
14704		1,000			
14027		5,000			
18365		50,000			
21985	Toyopearl Q-600C AR	100	50-150	0.15-0.20	>120
21986		250			
21987		1,000			
21988		5,000			
21989		50,000			
21854	Toyopearl GigaCap Q-650M	100	50-100	0.10-0.20	>170
21855		250			
21856		1,000			
21857		5,000			
21858		50,000			
19823	Toyopearl SuperQ-650S	25	20-50	0.20-0.30	105-155
17223		250			
17224		1,000			
17225		5,000			
19679		50,000			
43205	Toyopearl SuperQ-650M	100	40-90	0.20-0.30	105-155
17227		250			
17228		1,000			
17229		5,000			
21311		50,000			
43275	Toyopearl SuperQ-650C	100	50-150	0.20-0.30	105-155
17231		250			
17232		1,000			
17233		5,000			
19804	Toyopearl DEAE-650S	25	20-50	0.08-0.12	25-35
07472		250			
14692		1,000			
07973		5,000			
21483		50,000			



ION EXCHANGE CHROMATOGRAPHY

► ORDERING INFORMATION

ANION EXCHANGE RESINS:

TOYOPEARL BULK MEDIA

PART #	PRODUCT DESCRIPTION	CONTAINER SIZE (mL)	PARTICLE SIZE (µm)	ION EXCHANGE CAPACITY (meq/mL RESIN)	TYPICAL CAPACITY (mg BSA/mL RESIN)
43201	Toyopearl DEAE-650M	100	40-90	0.08-0.12	25-35
07473		250			
14693		1,000			
07974		5,000			
18367		50,000			
07988	Toyopearl DEAE-650C	250	50-150	0.05-0.11	25-35
14694		1,000			
07989		5,000			

TSKgel BULK RESIN

43383	TSKgel SuperQ-5PW (20)	25	15-25	0.12-0.18	52-88
18535		250			
18546		1,000			
18547		5,000			
43283	TSKgel SuperQ-5PW (30)	25	20-40	0.12-0.18	52-88
18536		250			
18548		1,000			
18549		5,000			
43381	TSKgel DEAE-5PW (20)	25	15-25	0.05-0.11	25-45
14710		250			
14711		1,000			
18436		5,000			
43281	TSKgel DEAE-5PW (30)	25	20-40	0.05-0.11	20-40
14712		250			
14713		1,000			
18370		5,000			

CATION EXCHANGE RESINS:

TOYOPEARL BULK MEDIA

PART #	PRODUCT DESCRIPTION	CONTAINER SIZE (mL)	PARTICLE SIZE (µm)	ION EXCHANGE CAPACITY (meq/mL RESIN)	TYPICAL CAPACITY (mg LYSOZYME/mL RESIN)
21833	Toyopearl GigaCap S-650M	100	50-100	0.14-0.18	80-120
21834		250			
21835		1,000			
21836		5,000			
21837		50,000			
21946	Toyopearl GigaCap CM-650M	100	50-100	0.17-0.28	>110 (γ-globulin)
21947		250			
21948		1,000			
21949		5,000			
21950		50,000			
43272	Toyopearl SP-550C	100	50-150	0.14-0.18	80-120
14028		250			
14705		1,000			
14029		5,000			
18366		50,000			

ION EXCHANGE CHROMATOGRAPHY



ORDERING INFORMATION


CATION EXCHANGE RESINS:
TOYOPEARL BULK MEDIA

PART #	PRODUCT DESCRIPTION	CONTAINER SIZE (mL)	PARTICLE SIZE (µm)	ION EXCHANGE CAPACITY (meq/mL RESIN)	TYPICAL CAPACITY (mg LYSOZYME/mL RESIN)
19822	Toyopearl SP-650S	25	20-50	0.13-0.17	40-60
08437		250			
14698		1,000			
08438		5,000			
21477		50,000			
43202	Toyopearl SP-650M	100	40-90	0.13-0.17	40-60
07997		250			
14699		1,000			
07998		5,000			
18369		50,000			
07994	Toyopearl SP-650C	250	50-150	0.12-0.18	35-55
14700		1,000			
07995		5,000			
19803	Toyopearl CM-650S	25	20-50	0.08-0.12	30-50
07474		250			
14695		1,000			
07971		5,000			
43203	Toyopearl CM-650M	100	40-90	0.08-0.12	30-50
07475		250			
14696		1,000			
07972		5,000			
19839		50,000			
07991	Toyopearl CM-650C	250	50-150	0.05-0.11	25-45
14697		1,000			
07992		5,000			
19329		50,000			
21804	Toyopearl MegaCap II SP-550EC	100	100-300	0.14-0.18	60-90*
21805		250			
21806		1,000			
21807		5,000			
21808		50,000			

TSKgel BULK MEDIA

PART #	PRODUCT DESCRIPTION	CONTAINER SIZE (mL)	PARTICLE SIZE (µm)	ION EXCHANGE CAPACITY (meq/mL RESIN)	TYPICAL CAPACITY (mg INSULIN/mL RESIN)
21976	TSKgel SP-3PW (30) NEW	25	20-40	0.07- 0.22	≥ 65
21977		250			
21978		1,000			
21979		5,000			
43382	TSKgel SP-5PW (20)	25	15-25	0.06-0.12	20-40
14714		250			
14715		1,000			
18435		5,000			
43282	TSKgel SP-5PW (30)	25	20-40	0.06-0.12	20-40
14716		250			
14717		1,000			
18384		5,000			

* Adsorption capacity for insulin: 90-120 mg/mL resin



HIC HYDROPHOBIC INTERACTION CHROMATOGRAPHY

HIC PRODUCTS

- Toyopearl SuperButyl-550
- Toyopearl Hexyl-650
- Toyopearl Butyl-650
- Toyopearl Phenyl-650
- Toyopearl Butyl-600
- Toyopearl Phenyl-600
- Toyopearl PPG-600
- Toyopearl Ether-650


- TSKgel Phenyl-5PW
- TSKgel Ether-5PW

≡ **TOSOH FACT**

Tosoh Bioscience offers a range of technical support services to our TSKgel, ToyoScreen, and Toyopearl chromatography products.

Whether you need help developing an HPLC assay for the analysis of a new therapeutic target, want to know how to monitor drug metabolites in the human body or need regulatory files to support a submission to the FDA, our technical support specialists will provide assistance in all of these areas and more.

We offer on-site training and application-specific seminars and are committed to providing prompt and courteous service for these and other requests.





HYDROPHOBIC INTERACTION CHROMATOGRAPHY

TOYOPEARL RESINS FOR HYDROPHOBIC INTERACTION CHROMATOGRAPHY

Hydrophobic interaction chromatography (HIC) is a powerful tool for the process purification of biomolecules. The technique utilizes the accessible hydrophobic regions located on protein surfaces and their interactions with a weakly hydrophobic stationary phase. HIC is an excellent complement to ion exchange (IEC) and size exclusion chromatography (SEC) particularly when protein isoforms exist or when feedstock impurities are of similar isoelectric point or molecular weight. The selectivity differences exploited by HIC can also be used after affinity separations in which closely related proteins with similar recognition sites are not distinguishable by the affinity ligand.


HOW DOES HIC WORK?

Proteins and other molecules with hydrophobic surfaces are attracted to the hydrophobic ligands of both reversed phase (RPC) and HIC resins. RPC resins have higher surface coverage and/or more hydrophobic ligand compared to HIC resins. Because of this, in a RPC separation the target binding readily occurs in an aqueous solution, and desorption is promoted by the addition of an increasing amount of organic solvent.

In HIC, proteins are bound to the resin by employing an aqueous high salt mobile phase. The salt conditions contribute to a lyotropic effect which allows the proteins to bind to the lower surface coverage of a hydrophobic ligand. Proteins are eluted by the simple technique of decreasing the salt concentration. Most therapeutic targets are eluted in a low salt or a no salt buffer.

During elution the energy of interaction for a HIC step is less than that of a RP step. One means of gauging the relative binding energy between the two techniques is to measure the surface tension of the two sets of binding and elution conditions. Figure 1 provides a comparison of the surface tension generated by HIC and RPC elution systems. Since HIC separates under milder eluting conditions, biological activity is typically retained.

FIGURE 1 SURFACE TENSION OF AQUEOUS SOLUTIONS USED IN HIC & RPC



Mode	Gradient (typical)	Δ Surface tension (erg/cm ²)
HIC	1.8 to 0 mol/L (NH ₄) ₂ SO ₄ / aqueous buffer	4
RPC	10 to 50% ACN/ 0.1%TFA	23

C. Horvath et. al., Separation Processes in Biotechnology, (J. Asenjo, Ed.) 9, 447 (1990) Marcel Dekker

FEATURES

➤ hydrophilic polymer resin matrix

➤ good mechanical stability

BENEFITS

- robust chemical stability between pH 1 - 13
- temperature range 4 - 60°C
- autoclavable at 121°C
- compatible with organic solvents
- constant bed volume over a wide range of salt concentrations
- low non specific protein binding
- superior protein recovery
- excellent flow characteristics in large industrial size columns
- direct scale-up from TSKgel HIC HPLC columns

HYDROPHOBIC INTERACTION CHROMATOGRAPHY



FIVE DIFFERENT HYDROPHOBIC SURFACES AND SELECTIVITIES

Tosoh Bioscience offers five HIC ligands featuring different degrees of hydrophobicity and selectivity. The hydrophobicity of Toyopearl HIC resins increases through the ligand series: Ether, PPG (polypropyleneglycol), Phenyl, Butyl, and Hexyl. (Figure 2)

Coordinating the hydrophobicity of the therapeutic target to the resin hydrophobicity is critical for the best overall purification performance. Too hydrophobic a resin for a given protein can result in its irreversible binding to the resin or a loss of enzymatic activity. Table I and II show typical mass recovery and biological activity recovery data for Toyopearl HIC resins.

An optimum HIC process step will balance high dynamic binding capacity, adequate selectivity, good mass recovery and retention of biological activity. The wide range of Toyopearl selectivities enables a developer to optimize protein separations at the extremes of the hydrophobic spectrum. Highly retentive Toyopearl Hexyl-type and Toyopearl Butyl-type resins are used to separate hydrophilic proteins. These two resins should also be considered for separations requiring a low salt environment. Toyopearl Ether-type resin is used for the purification of very hydrophobic targets such as certain monoclonal antibodies and membrane proteins. These proteins may bind irreversibly to other more hydrophobic resins. Toyopearl PPG-type and Toyopearl Phenyl-type phases complement the other HIC ligands available in the Toyopearl series and offer alternatives for mid-range hydrophobic proteins.

FIGURE 2

HIC LIGAND CANDIDATES




TABLE I

HIGH MASS RECOVERY (%) OF PROTEINS

	Toyopearl HIC resin		
	Ether-650M	Phenyl-650M	Butyl-650M
Bovine serum albumin	84	62	76*
α -chymotrypsinogen	96	88*	90
Cytochrome C	-	81*	87*
IgG	91	-	-
α -Lactalbumin	90	-	-
Lysozyme	94	92	85
Ovalbumin	83	88	73
Ribonuclease A	-	72*	82*

*Procedure: A 200 mL sample containing 200 mg of protein was loaded onto a 7.5 mm column and eluted with a 60 minute gradient of 1.8 mol/L (*1.5 mol/L) to 0.0 mol/L ammonium sulfate in 0.1 mol/L sodium phosphate (pH 7.0). The mass recovery was determined spectrophotometrically at UV 280 nm and 25°C.*

TABLE II

RECOVERY OF ENZYMATIC ACTIVITY OF PROTEINS

Toyopearl HIC resin	Protein	% Activity recovery
Phenyl-650	Phytochrome	79
Butyl-650	Halophilic protease	85
Butyl-650	Poly (3-hydroxybutyrate) depolymerase	88
Butyl-650	Aculeacin-A acylase	82
Butyl-650	Opine dehydrogenase	81



HYDROPHOBIC INTERACTION CHROMATOGRAPHY

The retention and selectivity of protein standards on Toyopearl HIC resins using the ToyoScreen process development columns are shown in Figure 3.

INFLUENCE OF SALT TYPE

In addition to the hydrophobicity of the ligand, the selectivity in HIC is influenced by the eluent salt type. Figure 4 demonstrates the effect of salt type on the resolution factor of different protein pairs.

The Hofmeister lyotropic salt series shown in Figure 5 ranks anions and cations by their ability to promote protein precipitation. Ions on the left are referred to as "lyotropic" while the ions on the right are called "chaotropic".

Lyotropic salts will precipitate or "salt out" proteins at high salt concentrations due to increased hydrophobic interaction, while chaotropic salts will promote protein denaturation at high salt concentrations. Figure 5 indicates that different salt systems may generate a variety of adsorption and desorption selectivities for each resin. This feature of HIC provides an additional parameter for the optimization of a process step.

FIGURE 5

HOFMEISTER LYOTROPIC SALT SERIES

For anions

$\text{SO}_4^{2-} > \text{HPO}_4^{2-} > \text{CH}_3\text{COO}^- > \text{halide} > \text{NO}_3^- > \text{ClO}_4^- > \text{SCN}^-$

For cations

$(\text{CH}_3)_4\text{N}^+ > \text{NH}_4^+ > \text{K}^+ > \text{Na}^+ > \text{Cs}^+ > \text{Li}^+ > \text{Mg}^{2+} > \text{Ca}^{2+} > \text{Ba}^{2+}$

Ammonium sulfate and sodium sulfate are the most commonly used salts in HIC. NaCl is often used as well.


FIGURE 6

METHACRYLIC BASE BEADS AVAILABLE FOR HIC

Pore size (Å)	50	125	400-500	750	1000	>1000	>1700
Product name							
Toyopearl HW-type	40	50	55	60	65	75	80
TSKgel PW-type	G1000	G2000	G4000	-	G5000	G6000	-

FIGURE 3

SCREENING OF TOYOPEARL HIC RESINS - STANDARD PROTEINS




Column: ToyoScreen (1 mL)

Mobile phase A: 0.1 mol/L phosphate buffer + 1.8 mol/L sodium sulfate (pH 7.0); Eluent B: 0.1 mol/L phosphate buffer (pH 7.0); Gradient: 30 min linear gradient from A to B; Flow rate: 1 mL/min; Inj. vol.: 50 µL; Detector: UV @ 280 nm; Samples: Ribonuclease A, Lysozyme, α-Chymotrypsinogen, 1 mg/mL

FIGURE 4

INFLUENCE OF SALT-TYPE ON RESOLUTION



Chromatography on a Toyopearl Butyl-substituted support
Column dimensions: 4.1 mm ID x 4 cm L

Mobile phase: Linear gradient, 20 min, 1.0 mol/L to 0 mol/L of indicated salt in 20 mmol/L phosphate buffer (pH 7.0);

Flow rate, 1 mL/min; Detector: UV @ 280 nm

J. Fausnaugh, L. Kennedy and F. Regnier, J. Chromatography 317, 141 (1984)

HYDROPHOBIC INTERACTION CHROMATOGRAPHY



PARTICLE SIZE OPTIMIZATION

Toyopearl and TSKgel PW-type methacrylic base beads incorporate the same polymer chemistry (see Figure 6 or refer to SEC section for a more detailed discussion) and are available in a variety of particle sizes:

100 μm	Toyopearl C-grade	Capture
65 μm	Toyopearl M-grade	Intermediate Purification
35 μm	Toyopearl S-grade	Intermediate Purification/Polishing
30 & 20 μm	TSKgel PW-type	High Resolution

In many cases analytical columns are available packed with 10 micron particles having the same selectivity as the process media. Figure 7 shows the variety of ligands and particle sizes available for HIC process-scale applications. This figure also depicts where each particle size is used in a chromatographic manufacturing train (i.e. for capture, intermediate purification, or polishing).

MECHANICAL STABILITY

The semi-rigid polymeric backbone of Toyopearl and TSKgel HIC resins permits high flow rates for maximum throughput and productivity. Toyopearl resins may be operated up to 3 bar and TSKgel PW-type resins may be operated up to 20 bar. The pressure-flow characteristics for each particle size grade of Toyopearl Phenyl-650 resins are shown in Figure 8. The beads are compatible with both fixed bed and dynamic axial compression columns.

CHEMICAL STABILITY

The polymeric structure of these products also makes them resistant to a wide range of pH conditions and ionic strengths. In addition, the hydroxylated surface of the base bead reduces non-specific binding of proteins.

FIGURE 7

Process step	Bead size	Process media
Capture	100 μm	Toyopearl SuperButyl-650C Toyopearl Hexyl-650C Toyopearl Butyl-650C Toyopearl Phenyl-650C
	65 μm	Toyopearl Butyl-600M Toyopearl Phenyl-600M Toyopearl PPG-600M
Intermediate Purification	35 μm	Toyopearl Butyl-650M Toyopearl Phenyl-650M Toyopearl Ether-650M
	30 μm	TSKgel Phenyl-5PW (30) TSKgel Ether-5PW (30)
Polishing	20 μm	TSKgel Phenyl-5PW (20) TSKgel Ether-5PW (20)
	10 μm	TSKgel Phenyl-5PW 2 mm ID x 7.5 cm L TSKgel Ether-5PW 2 mm ID x 7.5 cm L

Same selectivity HPLC columns are available for most process media



HYDROPHOBIC INTERACTION CHROMATOGRAPHY


PORE SIZE OPTIMIZATION

Most Toyopearl HIC products are derived from the versatile size exclusion resin, Toyopearl HW-65 (1,000 Å mean pore size), as the base bead for the majority of protein separations. But the pore size and accessible surface area of Toyopearl resins can be optimized for a given protein. More accessible surface area increases the dynamic binding capacity (DBC) of the bead for a particular therapeutic target. This has led to the development of two specialty lines of HIC materials with higher dynamic binding capacities.

HIGH DYNAMIC BINDING CAPACITY TOYOPEARL PROCESS RESINS


For monoclonal antibodies a pore size of 750 Å is optimum. A specially made base resin, Toyopearl HW-60, has this pore size. Three ligands are available on Toyopearl HW-60: polypropylene glycol (PPG), phenyl, and butyl. A comparison of their DBCs with Toyopearl Phenyl-650M resin is shown in Figure 9. The selectivities of Toyopearl Butyl-600M, Toyopearl PPG-600M and the Toyopearl Phenyl-600M resins, are shown in Figure 10.

FIGURE 9 BREAKTHROUGH CURVES OF POLYCLONAL IgG ON VARIOUS HIC RESINS




Column: 7.8 mm ID X 20 cm L; Sample: polyclonal human IgG
Binding buffer: 1 g/L IgG in 0.8 mol/L $(\text{NH}_4)_2\text{SO}_4$ + 0.1 mol/L sodium phosphate (pH 7.0); Linear velocity: 300 cm/hr
Temperature: 25 °C; Detector: UV @ 280 nm
DBC was calculated at 10% of breakthrough.

FIGURE 8 PRESSURE-FLOW CURVE FOR TOYOPEARL PHENYL-650 RESINS OF VARIOUS PARTICLE SIZE



Column: Toyopearl Phenyl-650C, M and S, 25 mm ID x 25 cm L
Mobile phase: 2 mol/L $(\text{NH}_4)_2\text{SO}_4$

FIGURE 10 COMPARISON OF TOYOPEARL 600M SERIES RESINS



Column: 7.5 cm ID X 7.5 cm L; Sample: 1 g/L RNase A (1), lysozyme (2) and α -chymotrypsinogen A (3)
Sample load: 100 μ L; Gradient: 60 min linear gradient from buffer A to B; Buffer A: 1.8 mol/L $(\text{NH}_4)_2\text{SO}_4$ + 0.1 mol/L sodium phosphate (pH 7.0); Buffer B: 0.1 mol/L sodium phosphate (pH 7.0)
Linear velocity: 136 cm/hr; Temperature: 25 °C
Detector: UV @ 280 nm

HYDROPHOBIC INTERACTION CHROMATOGRAPHY



Figure 11 compares the selectivities of the Toyopearl Phenyl-600M, Toyopearl Phenyl-650M and an Agarose Phenyl resin.

For smaller molecules such as lysozyme (12,000 Da) the even narrower pore diameter Toyopearl SuperButyl-550C resin (derived from the 500 Å pore diameter Toyopearl HW-55) is recommended. A comparison of the DBC of Toyopearl SuperButyl-550C with other Toyopearl HIC resins is shown in Figure 12.


The Toyopearl Phenyl-600M resin also has a high DBC for lysozyme (Figure 13). The engineered higher dynamic binding capacity of the 600 and 550 series HIC products for their specific targets and the selectivity differences induced by the smaller mean pore size of the respective beads can have a dramatic impact on process economics.

TOYOSCREEN PREPACKED COLUMNS FOR PROCESS DEVELOPMENT

ToyoScreen columns packed with the full range of our Toyopearl HIC products are available in 1 mL and 5 mL resin volumes. They provide a convenient way to screen different resins for both target retention and recovery. Multiple columns can be connected in series for additional separation. Please see the ordering information at the end of this section or contact us for more information on these products.

FIGURE 11


SELECTIVITY COMPARISON OF PHENYL-TYPE RESINS



Column: 7.5 cm ID X 7.5 cm L
 Sample: 1 g/L RNase A (1), lysozyme (2) and α-chymotrypsinogen A (3); Sample load: 100 µL; Gradient: 60 min. linear gradient from buffer A to B; Buffer A: 1.8 mol/L (NH₄)₂SO₄ + 0.1 mol/L sodium phosphate (pH 7.0); Buffer B: 0.1 mol/L sodium phosphate (pH 7.0)
 Linear velocity: 136 cm/hr; Temperature: 25°C; Detector: UV @ 280 nm

FIGURE 12

TYPICAL DYNAMIC BINDING CAPACITIES FOR LYSOZYME



Resin	Binding capacity (mg/mL) (10% Breakthrough)
Ether-650M	12.5
Phenyl-650M	27.5
Butyl-650M	32.2
Hexyl-650C	33.2
PPG-600M	44.2
SuperButyl-550C	54.3

Column size: 7.8 mm ID x 20 cm L; Sample: 1 mg/mL Lysozyme in 0.1 mol/L phosphate buffer + 1.8 mol/L sodium sulfate (pH 7.0); Linear Velocity: 100 cm/h; Detection: UV @ 280 nm



HYDROPHOBIC INTERACTION CHROMATOGRAPHY


MONOCLONAL ANTIBODIES

Hydrophobic interaction is a very useful technique for the purification of monoclonal antibodies. The diverse hydrophobic nature of mAbs is seen in Figure 14. This figure measures the hydrophobicity (using elution time as a surrogate measurement) of 51 different mouse IgGs on a TSKgel Phenyl-5PW analytical column. Some of the IgGs have elution times 2-3 times longer than others indicating greater hydrophobicity. The Toyopearl series of HIC ligands (Figure 2, page 33) with their different hydrophobicities gives chromatographic developers a range of options for finding the right ligand for their target molecule.

For a very hydrophobic mAb, such as mouse anti-chicken 14 kDa lectin, the less hydrophobic Toyopearl Ether ligand works quite well. The purification from ascites fluid (Figure 15) was performed with a 10 μm TSKgel Ether-5PW semi-preparative column. Identical selectivity for scale-up was found with corresponding 65 μm Toyopearl Ether-650M resin.

FIGURE 13

TOYOPEARL PHENYL-600M BREAKTHROUGH CURVE (LYSOZYME)



Binding capacity (mg/mL)
(10% Breakthrough)


Phenyl-600M
Phenyl-650M

58
27

Column: 7.8 mm ID x 20 cm L; Sample: 1 mg/mL lysozyme in 0.1 mol/L phosphate buffer (pH 7.0) + 1.8 mol/L $(\text{NH}_4)_2\text{SO}_4$
Linear velocity: 300 cm/h; Detector: UV @ 280 nm

FIGURE 14


HYDROPHOBIC DIVERSITY OF MOUSE MONOCLONALS



Plot of chromatographic elution times for 51 different mouse mAbs
Column: TSKgel Phenyl-5PW; Mobile phase: (A) 0.1 mol/L phosphate buffer containing 1.8 mol/L ammonium sulfate (pH 7.0); (B) 0.1 mol/L phosphate buffer (pH 7.0);
Flow rate: 1 mL/min; Gradient: (B) 0% (0 min)--0% (5 min)--100% (65 min) linear; Detector: UV @ 280 nm; Samples: 51 kinds of mouse monoclonal antibodies

FIGURE 15

PURIFICATION OF MABS FROM ASCITES FLUID



Column: A. TSKgel Ether-5PW, 7.5 mm ID x 7.5 cm L

B. Toyopearl Ether-650M, 7.5 mm ID x 7.5 cm L

Sample: anti-chicken 14 kDa lectin, diluted ascites fluid,
A. 1.5 mg in 100 μL ; B. 0.76 mg in 50 μL

Mobile phase: 60 min linear gradient from 1.5 mol/L to 0 mol/L $(\text{NH}_4)_2\text{SO}_4$ in 0.1 mol/L phosphate buffer (pH 7.0)

Linear velocity: 136 cm/h; Detection: UV @ 280 nm

HYDROPHOBIC INTERACTION CHROMATOGRAPHY



PROTEIN AGGREGATE REMOVAL

The larger pore Toyopearl products such as Toyopearl Butyl-650 and Toyopearl Phenyl-650 resins are very useful for protein aggregate separation and removal. Early work by Karger et al¹ in 1989 involving proteins and aggregates larger than 200 kDa demonstrated the effectiveness of HIC for this application.

GLYCOPROTEINS

Toyopearl HIC resins can purify glycoproteins, which often bind irreversibly to saccharide-based chromatographic media. Figure 16 shows the purification of a large glycoprotein on Toyopearl Butyl-650S resin.

DNA PLASMID PURIFICATION AND ENDOTOXIN REMOVAL

Toyopearl Hexyl-650C resin was successfully used for plasmid DNA purification by Cambrex, Baltimore, MD (US patent 6,953,686). Hexyl-650C was shown to be the most effective among HIC resins for endotoxin removal with capacities exceeding 2 million EU/mL of resin. Additionally, RNA and protein impurities were effectively eliminated. Hexyl-650C was also effective in separating the supercoiled and open circular forms of plasmid DNA.


REMOVAL OF MISFOLDED PROTEINS

Because misfolded proteins will generally be more hydrophobic than the native protein, Toyopearl Butyl-650M resin is used frequently for the removal of misfolded proteins. In many cases flow through chromatography can be accomplished under eluent conditions binding the misfolds while allowing the target to flow through the column.

Go to www.separations.eu.tosohbioscience.com, to check our database on the website for additional applications.

¹ Grinberg, N. Blanco, R. Yarmush, D. Karger, B. L. Protein Aggregation in High Performance Liquid Chromatography: Hydrophobic Interaction Chromatography of β -Lactoglobulin, Anal. Chem. 1989, 61, 514-520.

FIGURE 16 LARGE GLYCOPROTEIN PURIFIED ON TOYOPEARL BUTYL-650S



Column: Toyopearl Butyl-650S, 22 mm ID x 26 cm L
 Sample: crude protein from sea hare *Aplysia kurodai*;
 Mobile phase: multi-step $(\text{NH}_4)_2\text{SO}_4$ in 50 mmol/L phosphate buffer (pH 7.0); A. load & wash: 40% saturated $(\text{NH}_4)_2\text{SO}_4$; B. 20% saturated $(\text{NH}_4)_2\text{SO}_4$; C. 0% saturated $(\text{NH}_4)_2\text{SO}_4$; Detector: UV @ 280 nm



HYDROPHOBIC INTERACTION CHROMATOGRAPHY

ORDERING INFORMATION

TOYOSCREEN PROCESS DEVELOPMENT COLUMNS FOR HIC

PART #	PRODUCT DESCRIPTION	PACKAGE
21372	ToyoScreen Ether-650M	1 mL x 6 each
21373	ToyoScreen Ether-650M	5 mL x 6 each
21374	ToyoScreen Phenyl-650M	1 mL x 6 each
21375	ToyoScreen Phenyl-650M	5 mL x 6 each
21376	ToyoScreen Butyl-650M	1 mL x 6 each
21377	ToyoScreen Butyl-650M	5 mL x 6 each
21378	ToyoScreen Hexyl-650C	1 mL x 6 each
21379	ToyoScreen Hexyl-650C	5 mL x 6 each
21380	ToyoScreen PPG-600M	1 mL x 6 each
21381	ToyoScreen PPG-600M	5 mL x 6 each
21892	ToyoScreen Phenyl-600M	1 mL x 6 each
21893	ToyoScreen Phenyl-600M	5 mL x 6 each
21382	ToyoScreen SuperButyl-550C	1 mL x 6 each
21383	ToyoScreen SuperButyl-550C	5 mL x 6 each
21494	ToyoScreen Butyl-600M	1 mL x 6 each
21495	ToyoScreen Butyl-600M	5 mL x 6 each
21398	ToyoScreen HIC Mix Pack, (PPG-600M, Butyl-600M, Phenyl-650M, Butyl-650M, Phenyl-600M, Hexyl-650C)	1 mL x 6 Grades x 1 each
21399	ToyoScreen HIC Mix Pack, (PPG-600M, Butyl-600M, Phenyl-650M, Butyl-650M, Phenyl-600M, Hexyl-650C)	5 mL x 6 Grades x 1 each

TOYOSCREEN COLUMN ACCESSORIES

PART #	PRODUCT DESCRIPTION
21400	ToyoScreen Column Holder

TSKgel LABPAK

PART #	PRODUCT DESCRIPTION	CONTAINER SIZE (mL)	PARTICLE SIZE (µm)
43278	HICPAK PW (20) (Ether-5PW, Phenyl-5PW)	2 x 25 mL	10-30
43175	HICPAK PW (30) (Ether-5PW, Phenyl-5PW)	2 x 25 mL	20-40

TOYOPEARL HIC RESINS

PART #	PRODUCT DESCRIPTION	CONTAINER SIZE (mL)	PARTICLE SIZE (µm)	TYPICAL CAPACITY (mg LYSOZYME/mL RESIN)
43151	Toyopearl Ether-650S	25	20-50	10-30
16172		100		
16174		1,000		
16176		5,000		
19805	Toyopearl Ether-650M	25	40-90	10-30
16173		100		
16175		1,000		
16177		5,000		
21301	Toyopearl PPG-600M	25	40-90	45-55
21302		100		
21303		1,000		
21304		5,000		
21305		50,000		
21887	Toyopearl Phenyl-600M	25	40-90	45-65
21888		100		
21889		1,000		
21890		5,000		
21891		50,000		
43152	Toyopearl Phenyl-650S	25	20-50	30-50
14477		100		
14784		1,000		
14935		5,000		

HYDROPHOBIC INTERACTION CHROMATOGRAPHY





ORDERING INFORMATION

PART #	PRODUCT DESCRIPTION	CONTAINER SIZE (mL)	PARTICLE SIZE (μm)	TYPICAL CAPACITY (mg LYSOZYME/mL RESIN)
19818	Toyopearl Phenyl-650M	25	40-90	30-50
14478		100		
14783		1,000		
14943		5,000		
18364		50,000		
43126	Toyopearl Phenyl-650C	25	50-150	30-50
14479		100		
14785		1,000		
14944		5,000		
43153	Toyopearl Butyl-650S	25	20-50	30-50
07476		100		
14701		1,000		
07975		5,000		
18826		50,000		
19802		Toyopearl Butyl-650M		
07477	100			
14702	1,000			
07976	5,000			
18355	50,000			
43127	Toyopearl Butyl-650C	25	50-150	30-50
07478		100		
14703		1,000		
07977		5,000		
21448	Toyopearl Butyl-600M	25	40-90	40-60(g/L (γ-globulin))
21449		100		
21450		1,000		
21451		5,000		
21452		50,000		
19955	Toyopearl SuperButyl-550C	25	50-150	52-70
19956		100		
19957		1,000		
19958		5,000		
19959		50,000		
44465	Toyopearl Hexyl-650C	25	50-150	30-50
19026		100		
19027		1,000		
19028		5,000		
TSKgel 5PW HIC RESINS FOR HIGH RESOLUTION				
43276	TSKgel Ether-5PW (20)	25	10-30	10-30
16052		250		
16053		1,000		
18437		5,000		
43176	TSKgel Ether-5PW (30)	25	20-40	10-30
16050		250		
16051		1,000		
18439		5,000		
43277	TSKgel Phenyl-5PW (20)	25	10-30	10-30
14718		250		
14719		1,000		
18438		5,000		
43177	TSKgel Phenyl-5PW (30)	25	20-40	10-30
14720		250		
14721		1,000		
17210		5,000		



TOSOH

AFFINITY
CHROMATOGRAPHY




AFC AFFINITY CHROMATOGRAPHY

AFC PRODUCTS

- Toyopearl Resins for Protein A Affinity Chromatography - NEW
 - Toyopearl Resins for Affinity Ligand Coupling
 - Activated
 - Reactive
 - Toyopearl Resins with Group Specific
 - ToyoScreen Process Development Columns

TOSOH FACT

Tosoh Bioscience GmbH offers workshops on chromatography in process development & production. These three day courses deliver a comprehensive background to bioprocess chromatographic purification as an integrated process technique. The workshops provide a balance of effective presentations and practical hands-on experience with process scale and methods development equipment under the guidance of qualified and knowledgeable Tosoh Bioscience technical experts. Every year more than 30 participants use this opportunity to broaden or refresh their chromatographic and downstream processing knowledge.





AFFINITY CHROMATOGRAPHY

TOYOPEARL AFFINITY RESINS FOR PROCESS SCALE CHROMATOGRAPHY

There are many custom designed affinity ligands available to the chromatographer. Toyopearl affinity chromatography resins are functionalized with chemically active groups or group-specific ligands. Resins with activated functional groups are ready to directly couple a protein or other ligand. Resins with reactive groups require carbodiimide coupling or reductive amination to achieve a stable covalent linkage. The latest development in Toyopearl affinity resins is the high capacity Protein A affinity resin Toyopearl AF-rProtein A-650F developed to increase productivity in antibody purification.

PRESSURE-FLOW CHARACTERISTICS AND PHYSICAL/CHEMICAL STABILITY

Toyopearl resins remain dimensionally stable within wide extremes of pH and ionic strength. Moreover, the semi-rigid Toyopearl particles do not distort under flow rates that generate up to 3 bar pressure. These properties of the resins combined with the narrow particle size distributions result in superior pressure-flow characteristics for the packed Toyopearl bed. Linear velocities of 300 - 500 cm/h generate a pressure of between 1 and 2 bar in a 20 cm length bed. Achievement of high linear velocities at relatively low pressure enables high throughput production scale chromatography using equipment with moderate pressure limitations. Sanitization or cleaning may be conducted with up to 0.5 mol/L NaOH or 0.5 mol/L HCl depending upon the ligand. In affinity chromatography, in particular, the choice of cleaning agent will be largely dependent upon the chemical stability of the ligand, rather than that of the base resin.

➤ **TABLE I**

Protein A Affinity	Activated resin	Reactive resin	Group specific
Protein A	AF-Tresyl AF-Epoxy	AF-Amino AF-Carboxy AF-Formyl	AF-Blue HC AF-Red AF-Chelate AF-Heparin HC

➤ FEATURES

- Active, reactive and group specific resins
- New high capacity AF-rProtein A resin
- Standard 1.000 Å pore size
- Porous, hydrophilic polymer matrix
- High mechanical stability

TOYOPEARL RESIN FOR PROTEIN A AFFINITY CHROMATOGRAPHY


Toyopearl AF-rProtein A-650F is an affinity resin designed for efficient and robust large-scale purification of monoclonal antibodies. A newly developed recombinant protein A ligand (*E. coli*) was linked to the well proven methacrylic polymer backbone of Toyopearl media. The ligand is derived from one of the IgG-binding domains of the *staphylococcus aureus* protein A (Figure 1). Its amino acid sequence was optimized in order to increase its stability towards alkaline solutions.

Multipoint attachment of the ligand to the Toyopearl matrix further enhances chemical and thermal stability of the resin. In practice this pays off for a low level of protein A leaching and also for a high resistance to alkaline solutions employed in cleaning-in-place (CIP) procedures.

Toyopearl AF-rProtein-650F binds human and mouse immunoglobulin G, and Fab fragments.

➤ **FIGURE 1**

RECOMBINANT PROTEIN A DERIVED LIGAND



➤ BENEFITS

- Broad range of applications
- Efficient antibody purification
- High capacity for large biopolymers
- Suitable for laboratory and process scale purifications
- Constant bed volume over a wide range of buffer composition

AFFINITY CHROMATOGRAPHY



HIGH BINDING CAPACITY AT SHORT RESIDENCE TIME


The particle size of 45 µm and the ligand density of Toyopearl AF-rProtein A-650F were optimized in order to reach high dynamic binding capacities (DBC) for immunoglobulins. Typical static IgG binding capacity is > 45 mg/ml resin and typical dynamic IgG binding capacity at 10 % breakthrough is > 30 mg/mL resin at 2 minutes residence time (1 mg/mL protein load). Fast mass transfer kinetics support high binding capacities even when applying high titer feedstocks at high flow rates. IgG breakthrough curves (Figure 2) at various linear velocities demonstrate the high IgG DBC at high velocities and the superior kinetic performance of Toyopearl AF-rProtein A-650F.

HIGH CIP AND SANITIZATION STABILITY

The structure of the recombinant ligand and its multipoint attachment to the base matrix enhances the stability of Toyopearl AF-rProtein A-650F in 0.1-0.5 M NaOH. The dynamic binding capacity remains high after repeated CIP cycles. After more than 150 CIP cycles with 0.1 M NaOH at 16 min contact time per cycle more than 90 % of initial dynamic binding capacity was retained (Figure 3). When performing cleaning-in-place with 0.5 M NaOH about 80 percent of IgG binding capacity remain after 50 cycles with 17 min contact time.


Toyopearl AF-rProtein A-650F is also stable in ethanol, 6 M urea, 6 M guanidinium chloride, and 1% phosphoric acid, respectively. Static binding capacity of the resin is not impaired when heated for 30 minutes to temperatures of up to 90°C. It can be stored at room temperature at production site. Recommended conditions for long term storage are a storage solution of 20 % ethanol and temperature of 4 - 8 °C.

FIGURE 2
DYNAMIC BINDING CAPACITY



Breakthrough curves for h-IgG loading (polyclonal, 10 mg/ml)
 Typical DBC at 10% breakthrough: 30,5 mg/mL @ 100 cm/hr (3 min residence time) - 24 mg/mL @ 200 cm/hr (1.5 min residence time)
 Column: 5 mm ID x 5 cm L;
 Mobile phase: 20 mM sodium phosphate buffer pH 7.2 containing 150 mM NaCl; Sample conc.: 10 mg/mL; Residence time: 1.5, 2.0, 2.5, 3.0 min

FIGURE 3
CLEANING-IN-PLACE STUDY WITH 0.1 M NaOH



Column: 5 mm ID x 5 cm L
 10 column volumes binding buffer pH 7.4
 5 column volumes elution buffer pH 3.0
 3 column volumes binding buffer containing 0.1 M NaOH, 16 min contact time
 3 column volumes binding buffer pH 7.4



AFFINITY CHROMATOGRAPHY

TOYOPEARL RESINS FOR AFFINITY LIGAND COUPLING

Toyopearl offers a spectrum of carefully selected affinity resins primed with activated or reactive groups which can be used to covalently attach almost any custom ligand. The structures of Toyopearl activated and reactive ligands are shown in Figure 4.


In general, Toyopearl AF-Tresyl-650M and Toyopearl AF-Formyl-650M resin are recommended for coupling proteins, while Toyopearl AF-Epoxy-650M resin is suited for coupling lower molecular weight ligands. Toyopearl AF-Amino-650M and Toyopearl AF-Carboxy-650M resins may be used for both.

Toyopearl affinity resins may be used in combinatorial chemistry or for solid phase synthesis of peptides and oligonucleotides because of their excellent stability in a variety of organic solvents and under extremes of pH.

FIGURE 4


ACTIVATED TOYOPEARL AFFINITY RESINS

Toyopearl AF-Tresyl-650M



Ligand Density: 80 $\mu\text{mol/g}$ (dry)


Toyopearl AF-Epoxy-650M



Ligand Density: 800 $\mu\text{mol/g}$ (dry)


REACTIVE TOYOPEARL AFFINITY RESINS

Toyopearl AF-Formyl-650M




Ligand Density: 60 $\mu\text{eq/mL}$

Toyopearl AF-Amino-650M



Ligand Density: 100 $\mu\text{mol/mL}$

Toyopearl AF-Carboxy-650M



Ligand Density: 100 $\mu\text{eq/mL}$

TABLE II

REPRESENTATIVE COUPLING DENSITIES FOR ACTIVATED AND REACTIVE TOYOPEARL MEDIA

Toyoperl resin	AF-Tresyl-650M	AF-Formyl-650M	AF-Amino-650M	AF-Carboxy-650M
Protein coupled (mg/mL resin)				
Soybean trypsin inhibitor	16	3.5	5.8	15
Protein A	1.9	-	-	-
Concanavalin A	13	-	-	-
α 1-Antitrypsin	12.3	-	-	-
α -Chymotrypsin	12.5	-	-	-
Myoglobin	12.4	-	-	-
Ovalbumin	-	2.5	6.7	0.8
Bovine serum albumin	12.4	14	19.2	3.3
Human IgG	10.0	15	6.7	11.7
Cytochrome C	-	5.8	3.3	7.5
Lysozyme	60	20	5.8	17.5
Coupling agent	not required	NaCNBH ₃	NaCNBH ₃ or Cabodiimidamide	Carbodiimide
Optimal pH	7.0 - 9.0	6.9 - 9.0	4.5 - 6.0	4.5 - 6.0

AFFINITY CHROMATOGRAPHY



ACTIVATED RESINS – READY FOR DIRECT LIGAND ATTACHMENT

Toyopearl AF-Tresyl-650M activated resin is highly reactive toward amine and thiol groups. It is provided in dry form, ready for reaction in buffered solutions containing protein or other ligand. Coupling is accomplished in neutral to slightly alkaline (pH 7 - 8) solution (Figure 5).

Under such conditions, even proteins of limited stability may be successfully coupled. Coupling leads to the formation of a highly stable secondary amine or thio-ether linkage. The optimized tresyl-density (ca. 20 $\mu\text{mol/mL}$ hydrated resin) is sufficient to provide substantial protein binding while avoiding excessive multi-point attachment and consequent impairment of ligand affinity/activity. Representative data are presented in Table II.

Toyopearl AF-Epoxy-650M activated resin, also packaged in dry form, has a high density of epoxy- functionality (ca. 800 $\mu\text{mol/mL}$). Under appropriate reaction conditions, this may be used for immobilization of proteins or low molecular weight ligands. It is particularly useful when high densities of low molecular weight ligands must be attached (Figure 6). Glutathione and glycine have, for example, been coupled at densities greater than 100 $\mu\text{mol/mL}$ hydrated resin.

Toyopearl AF-Epoxy-650M resin is a highly versatile starting material for conversion to other chemically active functional groups required in special applications. This resin may be readily activated to hydrazide-bearing materials. This is particularly useful for immobilization of carbohydrates or glycoproteins. Using the reaction sequences described, special ligands may be introduced onto this dimensionally stable, macroporous support.

FIGURE 5 TOYOPEARL AF-TRESYL COUPLING PROCEDURE




FIGURE 6 TOYOPEARL AF-EPOXY COUPLING PROCEDURE




FIGURE 7 TOYOPEARL AF-FORMYL COUPLING PROCEDURE





FIGURE 8 TOYOPEARL AF-AMINO COUPLING PROCEDURE





AFFINITY CHROMATOGRAPHY

REACTIVE RESINS - REQUIRE ACTIVATION FOR LIGAND ATTACHMENT

Ligands may be coupled to Toyopearl AF-Formyl-650M (aldehyde-bearing) resin under mild conditions exclusively using primary amines. The ligand is bound to the resin by a stable secondary amine linkage (Figure 7). Representative coupling capacities are shown in Table II.

A wide variety of industrial enzymes have been immobilized on aldehyde-bearing supports. Typically, these supports have been synthesized by industrial users by partial oxidation of polysaccharide supports (e.g. cellulose and agarose) or partial hydrolysis of polyacetals. In contrast, Toyopearl AF-Formyl-650M resin is a ready-to-use aldehyde support formulated from a dimensionally stable, macro-porous matrix. Consistent aldehyde content and physical properties are assured from batch to batch.


Toyopearl AF-Amino-650M resin may be used to couple ligands using their carboxyl groups (peptide bond formation) or aldehyde groups (reductive amination) as shown in Figure 8. Aldehyde groups may be present in a carbohydrate or glycoprotein ligand or may be introduced into the ligand by mild, periodate oxidation.

The optimized functional group density of Toyopearl AF-Amino-650M (100 $\mu\text{mol/mL}$) is ideal for coupling of either proteins or low molecular weight ligands. For example, lactose was coupled by reductive alkylation to yield a ligand density of ca. 30 $\mu\text{mol/mL}$ resin. Coupling densities for various proteins are given in Table II.

Toyopearl AF-Carboxy-650M resin provides another useful and milder approach for coupling to amino groups of proteins or low molecular weight ligands. The carbodiimide mediated coupling reaction produces an amide bond between ligand and support (Figure 9). Representative coupling densities are given in Table II.

FIGURE 9

TOYOPEARL AF-CARBOXY COUPLING PROCEDURE



TOYOPEARL RESINS WITH GROUP SPECIFIC LIGANDS

The structures of Toyopearl group specific ligands are shown in Figure 10.

Toyopearl AF-Chelate-650M

This resin is derivatized with iminodiacetic acid (IDA) at a concentration of ca. 20 $\mu\text{mol/mL}$. In typical applications, selected metal ions, most often Ca^{2+} , Ni^{2+} , Zn^{2+} , Co^{2+} and Cu^{2+} are bound to the support by stable chelation. The resultant metal ion-bearing resin binds to histidine and free cysteine containing sequences of a peptide or protein. Immobilized metal ion affinity chromatography (IMAC) has been used for purification of recombinant human growth factor, tissue plasminogen activator, glycoporphins, and whole cells.


Toyopearl AF-Blue HC-650M

Functionalized with Cibachron Blue F3G-A, Toyopearl AF-Blue HC-650M resin has excellent capacity for proteins, particularly albumin (Figure 10). In addition, this high capacity resin has improved caustic stability, reduced dye ligand leakage, and superior pressure-flow characteristics relative to more traditional agarose materials (Figure 11).


FIGURE 10

GROUP-SPECIFIC TOYOPEARL AFFINITY RESIN


Toyopearl AF-Blue HC




Toyopearl AF-Red



Toyopearl AF-Chelate



Toyopearl AF-Heparin



AFFINITY CHROMATOGRAPHY



Toyopearl AF-Red-650ML


Toyoscreen and Toyopearl AF-Red-650ML resins are functionalized with Procion Red HE-3B, (also known as Reactive Red 120). This resin is useful for the purification of nucleotide dependent enzymes, lipoproteins, plasminogen, peptides, hormones and cytotoxins.

These two dye-ligand resins are useful in binding/purification of nucleotide-dependent enzymes, albumin, cell growth factors, interferons, transferases, cyclases, and polymerases. Typical binding capacities are shown in Table III.

TABLE III
REPRESENTATIVE BINDING CAPACITIES FOR TOYOPEARL DYE-LIGAND AFFINITY MEDIA

Protein (mg/mL resin)	AF-Blue HC-650M	AF-Red-650ML
Hexokinase	3	-
Bovine serum albumin	16	-
Human serum albumin	18+-2.5	3.5+-1
Lactate dehydrogenase	27	11

FIGURE 11
COMPARISON OF HUMAN SERUM ALBUMIN BINDING CAPACITIES AT VARIOUS pH'S OF AF-BLUE HC-650M AND AGAROSE (BLUE FUNCTIONALIZED AGAROSE)



Conditions

- A 1.0 mL of adsorbent was washed with 10 mL of equilibration buffer (pH 4.0 and 5.0; 0.1 mol/L sodium acetate buffer, pH 6.0 and 7.0; 0.1 mol/L sodium phosphate buffer).
- A 5.0 mL of 1 % solution of human albumin dissolved in each equilibration buffer was charged onto the column.
- After 10 min, unbound albumin was eluted and the column was washed with 10 mL of each equilibration buffer.
- Adsorbed human albumin was eluted with 0.1 mol/L sodium phosphate buffer at pH 7.0 containing 2.0mol/L sodium chloride (desorption buffer) and 10 mL fractions were collected.
- Human albumin content was measured spectrophotometrically by using E 0.1% at 280 nm = 0.55.

Toyopearl AF-Heparin HC-650M


Heparin is a linear and highly sulfated glycosaminoglycan which has anti-coagulant properties. Due to its polyanionic nature, heparin interacts with a wide range of biomolecules including plasma components, lipoprotein lipase, collagenase, and DNA polymerase.

Immobilized heparin is widely used as an adsorbent in affinity chromatography for the purification of biological substances. Toyopearl AF-Heparin HC-650M resin is a high capacity, affinity adsorbent with excellent chemical stability.

TOYOSCREEN PREPACKED COLUMNS FOR PROCESS DEVELOPMENT

ToyoScreen columns are available in 1 mL and 5 mL resin volumes of the following Toyopearl AFC products: Toyopearl AF-Chelate 650M, Toyopearl AF-Heparin HC-650M, Toyopearl AF-Blue HC-650M and Toyopearl AF-Red-650M. They provide a convenient way to perform early resin screening for both target retention and recovery. Multiple columns can be connected in series for additional separation. Please see the ordering information at the end of this section or contact us to request more information on our ToyoScreen offerings.

FIGURE 12
COMPARATIVE DYE LEAKAGE STUDY OF AF-BLUE HC-650M AND COMPETITOR BLUE @ 25°C (AFTER 24 HOURS)



Conditions: 200 mg of each material was suspended in 4 mL of solvent and incubated at 25°C and shaking for 24 h. The absorption at 620 nm of the supernatants were measured after appropriate adjustment to neutrality with known volume of acid or base. Dye concentrations were estimated assuming a molar extinction coefficient of 12,750 (L/M cm).



AFFINITY CHROMATOGRAPHY

ORDERING INFORMATION

TOYOSCREEN PROCESS DEVELOPMENT COLUMNS FOR AFC

PART #	PRODUCT DESCRIPTION	PACKAGE
22809	ToyoScreen AF-rProtein A-650F NEW	1 mL x 5 each
22810	ToyoScreen AF-rProtein A-650F NEW	5 mL x 1 each
22811	ToyoScreen AF-rProtein A-650F NEW	5 mL x 5 each
21384	ToyoScreen AF-Chelate-650M	1 mL x 6 each
21385	ToyoScreen AF-Chelate-650M	5 mL x 6 each
21386	ToyoScreen AF-Blue HC-650M	1 mL x 6 each
21387	ToyoScreen AF-Blue HC-650M	5 mL x 6 each
21388	ToyoScreen AF-Red-650ML	1 mL x 6 each
21389	ToyoScreen AF-Red-650ML	5 mL x 6 each
21390	ToyoScreen AF-Heparin HC-650M	1 mL x 6 each
21391	ToyoScreen AF-Heparin HC-650M	5 mL x 6 each

TOYOSCREEN COLUMN ACCESSORIES

PART #	PRODUCT DESCRIPTION
21400	ToyoScreen Column Holder

TOYOPEARL LABPAK

PART #	PRODUCT DESCRIPTION	CONTAINER SIZE (mL)	PARTICLE SIZE (µm)
43400	AFFIPAK ACT (AF-Epoxy-, AF-Tresyl-650M)	2 x 5 g*	65
43410	AFFIPAK (AF-Amino-, AF-Carboxy-, AF-Formyl-650M)	3 x 10 mL	65

TSKgel RESIN

PART #	PRODUCT DESCRIPTION	CONTAINER SIZE (mL)	PARTICLE SIZE (µm)
16208	Tresyl-5PW (10)	2 g	10

TOYOPEARL AFFINITY CHROMATOGRAPHY RESIN

GROUP SPECIFIC RESINS

PART #	PRODUCT DESCRIPTION	CONTAINER SIZE (mL)	TYPICAL CAPACITY
22803	Toyopearl AF-rProtein A-650F NEW	10	40 mg/mL (IgG)
22804		25	
22805		100	
22806		1,000	
22807		5,000	
22808		50,000	

PROTEIN A IMMUNOASSAYS & STANDARDS

PART #	PRODUCT DESCRIPTION
22815	Protein A-R28 ELISA Kit NEW
22836	Protein A-R28 STD 0.5 mL (10 mg/L) NEW

AFFINITY CHROMATOGRAPHY



ORDERING INFORMATION

GROUP SPECIFIC RESINS

PART #	PRODUCT DESCRIPTION	CONTAINER SIZE (mL)	TYPICAL LIGAND DENSITY	TYPICAL CAPACITY
19688	Toyopearl AF-Blue-HC-650M	25	15 µmol/mL	minimum 18 mg/mL
19689		100		
19690		1,000		
19691		5,000		
08651	Toyopearl AF-Red-650ML	25	5 µmol/mL	2.5 - 4.5 mg/mL (HSA)
19801		100		
42102		1,000		
14475	Toyopearl AF-Chelate-650M	25	25-45 µeq/mL	-
19800		100		
14907		1,000		
14908		5,000		
20030	Toyopearl AF-Heparin HC-650M	10	-	5 mg/mL (AT III)
20031		100		
20032		1000		
20033		5000		

TOYOPEARL AFFINITY CHROMATOGRAPHY RESIN

REACTIVE RESINS

PART #	PRODUCT DESCRIPTION	CONTAINER SIZE (mL)	TYPICAL LIGAND DENSITY	TYPICAL CAPACITY
43411	Toyopearl AF-Amino-650M	10	70-130 µeq/mL	-
08002		25		
08039		100		
18074		1,000		
18316		5,000		
43412	Toyopearl AF-Carboxy-650M	10	80-120 µeq/mL	-
08006		25		
08041		100		
18827		1,000		
18828		5,000		
43413	Toyopearl AF-Formyl-650M	10	40-70 µeq/mL	-
08004		25		
08040		100		
17396		1,000		
17397		5,000		

ACTIVATED RESINS


PART #	PRODUCT DESCRIPTION	CONTAINER SIZE	TYPICAL LIGAND DENSITY	TYPICAL CAPACITY
43402	Toyopearl AF-Epoxy-650M	5 g*	600 - 1000 µeq/g	-
08000		10 g*		
08038		100 g*		
18315		1,000 g*		
14471	Toyopearl AF-Tresyl-650M	5 g*	80 µmol/mL	-
14472		100 g*		
14906		1,000 g*		

Conditions: All Toyopearl affinity resins are provided at a particle size of 65 µm. This particle size is ideal for both small and large scale separations.

*1 g yields approximately 3.5 mL of hydrated resin.




TOSOH



CUSTOM DESIGNED RESINS

CUSTOM RESINS

CUSTOM DESIGNED RESINS



TOSOH FACT

Tosoh Bioscience provides solutions for today's biological purification needs. In fact, some of the first commercial HIC products were manufactured by Tosoh. We take pride in our ability to design new products based on existing chemistries to solve specific customer applications.

We encourage you to have a confidential discussion with us about your specific needs. Whether it is a surface modification of an existing product or the creation of a new one, we encourage you to call on us to meet your needs for a customized solution.



CUSTOM DESIGNED RESINS

OVERVIEW

Occasionally special purification needs require special adaptations to existing chromatographic resins for optimized operation of your process. This section explains in general terms Tosoh's custom resin development procedure.

TOSOH'S RESIN INNOVATION PROGRAM (TRIP)

Tosoh Corporation of Japan has introduced a new resin innovation program. Ideal candidates for this program are drug targets that have been selected for clinical trials which, if successful in the clinic, will be commercially manufactured. This program is not available for early stage laboratory research or the isolation of small amounts of new therapeutic proteins.

NEW PROJECT APPROVAL FOR TRIP

To receive a "Go" status for TRIP some evaluation is needed. During preliminary discussions, prior to the generation of scouting samples, an assessment is made of whether the project fits into Tosoh's scaffolding capabilities. If a fit is established, the projected liter volume for the new resin is reviewed. Decisions are based on the quantity projected for use after therapeutic approval and launch, but for high value added resins, projected commercial liter volumes of less than 100 liters have been approved.

NEW RESIN DEVELOPMENT TIME LINE

Tosoh Bioscience and Tosoh Corporation of Japan have the resources to develop custom resin products. The custom resin optimization timeline closely mirrors a customer's need to produce products for clinical trials. Only with detailed sample evaluation data from the customer, chemistry optimization and manufacturing lot release criteria can be established by Tosoh so a close cooperation between the customer and Tosoh is crucial.

TABLE 1

TRADITIONAL CHROMATOGRAPHIC LIGANDS

Anion Exchangers	Diethylaminoethyl (DEAE) Quaternary Anion Exchanger (QAE)
Cation Exchangers	Sulfopropyl (SP) Carboxymethyl (CM) Sulfonic Acid (S)
Hydrophobic Interaction	Hexyl Butyl Phenyl Polypropylene glycol (PPG) Ether

TABLE 2

TOSOH METHACRYLIC BASE BEADS USED FOR SEC

Pore size (Å)	50	125	400-500	750	1000	>1000	>1700
Product name							
Toyopearl HW-Type	40	50	55	60	65	75	80
TSKgel PW-Type	G1000	G2000	G4000		G5000	G6000	

TABLE 3

AVAILABLE PARTICLE AND PORE SIZE COMBINATIONS

Pore Diameter (angstroms)	Particle Size (microns)							
	200	100	75	65	35	30	20	15
>1700		★						
>1200	★	★						
1000		★		★	★	★	★	★
750		★	★	★				
400-500	★	★		★	★	★	★	★
125				★	★			
50	★	★		★	★			
	★ commercial	★ experimental						

7 x 8 = 56 possible combinations

CUSTOM DESIGNED RESINS



RESIN OPTIMIZATION (SEMI-CUSTOM COMMERCIAL RESINS)

In many cases the optimization of an already commercial resin can be accomplished by a more judicious selection of particle size, pore size, and currently available ligands (Table I). See the various combinations of these physical parameters in the Toyopearl SEC section of this catalog. When these products are combined with the comparable polymer chemistry TSKgel PW products (Table II), there are 49 different combinations (Table III) to act as a scaffold for an improved resin product. All kinds of modes can be optimized.

CUSTOM RESINS

In more difficult situations custom resins can be developed for a customer. Ligand screening, selection, density, length of spacer arm, degree of crosslinking and degree of grafting all add to the fundamental scaffold of the semi-custom resins. The objective in designing a certain resin is to develop a bead that has maximum accessible surface area with an optimized ligand density using an appropriate spacer arm which fits the conformational attachment requirements of the target. These parameters can have a dramatic effect on the resin's chromatographic and economic performance.

GETTING STARTED

To initiate the TRIP program contact your local Business Development Manager. A period of due diligence will be started where the appropriate information can be exchanged in confidence between us. Detailed communications on these topics requires an in effect two way Confidential Information Disclosure Agreement between Tosoh Bioscience GmbH and the customer. If the project looks promising to both parties several laboratory scouting samples will be prepared and sent to the customer site for evaluation.

COMMERCIALIZATION

Once a "Go" decision is made, the performance specifications are determined from the optimized scouting samples, the custom resin enters into Tosoh's "ISO 9001" product commercialization process. Several small manufacturing batches are then made to set lot release criteria. Regulatory support studies are initiated at this time.



YOU'VE SPECIALIST IN SEPARA RATION



MORE PRODUCTIVE MAB PURIFICATION

- EFFICIENT MAB CAPTURING FROM HIGH TITER FEEDSTOCK
- HIGH IgG CAPACITY OVER A WIDE RANGE OF VELOCITIES
- ROBUST, ALKALI-STABLE AFFINITY RESIN
- RIGID POLYMETHACRYLATE BASE MATRIX

TOYOPEARL AF-rPROTEIN A-650F IS A HIGH CAPACITY AFFINITY RESIN DESIGNED FOR EFFICIENT PURIFICATIONS OF MONOCLONAL ANTIBODIES. LOW LIGAND LEACHING AND HIGH MECHANICAL, THERMAL AND ALKALI STABILITY IMPROVE ROBUSTNESS AND PROCESS ECONOMICS.

TOSOH BIOSCIENCE

TOSOH BIOSCIENCE GMBH | ZETTACHRING 6 | 70567 STUTTGART, GERMANY
PHONE · + 49 (0) 711 · 13257-0 | FAX · + 49 (0) 711 · 13257-89
info.tbg@tosoh.com | www.tosohbioscience.com

FURTHER INFORMATION



For detailed Toyopearl packing instructions, request our **Toyopearl Instruction Manual**.

To get an overview about the whole range of TSKgel columns and small Toyopearl and TSKgel bulk media, please request our **Chromatography catalog**.

TOSOH

TOYOPEARL INSTRUCTION MANUAL

TOSOH BIOSCIENCE

TOSOH

Chromatography Catalog

YOUR SPECIALIST IN SEPARATION

TOSOH BIOSCIENCE

For a deeper insight into applications and all questions related to the practical use of TSKgel and Toyopearl, check out the website www.tosohbioscience.com and related catalogs or instruction manuals.

Our technical experts are happy to discuss your specific separation needs by phone: **+49 (0)711 13257-57** or techsupport.tb@tosoh.com



TOSOH

TOSOH BIOSCIENCE

Zettachring 6, 70567 Stuttgart, Germany
Tel: +49 (0)711 13257 0 Fax: +49 (0)711 13257 89
info.tb@tosoh.com www.tosohbioscience.com